



GC AND GC/MS



In this Chapter

252 GC and GC/MS Maintenance Schedule

Parts and Supplies

- 254** GC General Supplies
- 280** GC Instrument Parts and Supplies
- 318** GC/MS Parts and Supplies
- 337** GC Parts and Supplies for Varian Instruments

Agilent J&W GC and GC/MS Columns

- 359** Ultra Inert Columns
- 362** Guard Columns
- 364** Low-bleed GC/MS Columns
- 392** Polysiloxane Polymers Columns
- 430** Polyethylene Glycol (PEG) Columns
- 442** PLOT Columns
- 461** Special Application Columns
- 524** Packed GC Columns
- 531** Test Standards
- 532** Fused Silica Tubing
- 559** GC Column Troubleshooting and Maintenance

GC and GC/MS Applications

- 568** Environmental Applications
- 617** Foods, Flavors and Fragrance Applications
- 641** Industrial Chemical Applications
- 671** Life Science Applications
- 686** Petroleum Applications





GC and GC/MS Maintenance Schedule

Item	Typical Schedule	Actions/Comments
Gas Management		
Gas purifiers (carrier gas and detector gas)	Every 6 to 12 months	Replacement schedule is based on capacity and grade of gas. In general, replace non-indicating traps every 6 to 12 months or when indicating traps start to change color. Replace indicating traps when indicating material is starting to change color.
Internal split vent trap	Every 6 months*	Replace to prevent material backing up into EPC control and to avoid costly repair.
External split vent trap	Every 6 months*	Replace to prevent sample analytes from escaping into the laboratory environment.
Flow meter calibration	Every 1 to 2 years	Re-calibrate electronic flow meters – follow recommended schedule for the unit (shown on the calibration certificate).
Sample Introduction and Inlets		
Syringes and/or syringe needles	Every 3 months*	Replace syringe if dirt is noticeable in the syringe, if it cannot be cleaned, if the plunger doesn't slide easily, or if clogged. Replace needle if septa wear is abnormal or the needle becomes clogged.
Inlet liner	Weekly*	Check often. Replace when dirt is visible in the liner or if chromatography is degraded.
Liner O-rings	Monthly*	Replace with every liner change.
Inlet septum	Daily*	Check often. Replace when signs of deterioration are visible (gaping holes, fragments in inlet liner, poor chromatography, low column pressure, etc.)
Inlet hardware	Every 6 months Every year	Check for leaks and clean. Check parts and replace when parts are worn, scratched, or broken.
Inlet gold or stainless steel seal	Monthly*	For highest level of reproducibility, change inlet seal with every liner change, but minimally replace monthly or when scratched, corroded, or if there is build-up of non-volatile sample components.

*Schedule is an approximation of average usage requirements. Frequency may vary widely based upon application and sample type.

(Continued)

GC and GC/MS Maintenance Schedule		
Item	Typical Schedule	Actions/Comments
Columns		
Front-end maintenance	Weekly-monthly*	Remove 1/2 to 1 m from the front of the column when experiencing chromatographic problems (peak tailing, decreased sensitivity, retention time changes, etc.). Replace inlet liner and septum, and clean inlet as necessary. Guard column may be useful for increasing column lifetime.
Solvent rinse	As needed	Perform when chromatography degradation is due to column contamination. Only for bonded and cross-linked phases.
Replacement	As needed	Replace when trimming and/or solvent rinsing no longer restore chromatographic performance.
Ferrules	As needed	Replace when changing columns and inlet/detector parts.
Detectors		
FID/NPD jets and collector	As needed	Clean when deposits are present. Replace when they become scratched, bent, or damaged, or when having difficulty lighting FID or keeping flame lit.
NPD bead	As needed	Replace when signal drifts or there is a dramatic change in sensitivity.
FID	Every 6 months	Measure hydrogen, air, and makeup gas flows.
TCD	As needed	Thermally clean by "baking-out" when a wandering baseline, increased noise, or a change in response is present. Replace when thermal cleaning does not resolve the problem.
ECD	Every 6 months or as needed	Wipe test. Thermally clean by "baking-out" when baseline is noisy, or the output value is abnormally high. Replace when thermal cleaning does not resolve the problem.
FPD	Every 6 months or as needed	Measure hydrogen, air, and makeup gas flows. Clean/replace FPD windows and seals when detector sensitivity is reduced.
NCD and SCD	Every 3 months*	Change pump oil, oil coalescing filter and chemical trap.
Mass Selective Detectors		
Tune MSD	As needed	Keep plenty of PFTBS (P/N 05971-60571) on hand.
Check the calibration vial	Every 6 months	Vial can be refilled without venting the system.
Replace the foreline pump oil	Every 6 months	Check the fluid weekly. Change when the fluid becomes discolored or every 6 months.
Replace the diffusion pump fluid	Every year or as needed	Check the fluid weekly. Too little fluid will cause the pump to run at a higher temperature, resulting in degradation and loss of high vacuum. Change when the fluid is discolored or contains particulates.
Clean the ion source	As needed	Clean when performance deteriorates to remove contamination and to restore the electrostatic properties of the ion lens system. Replace scratched parts to maintain optimal performance.

*Schedule is an approximation of average usage requirements. Frequency may vary widely based upon application and sample type.

PARTS AND SUPPLIES

NEW!

Bulk GC Supplies



Single taper splitless liner, no wool, 5190-2270



Certified gold plated seal kit, 5190-2209



Certified non-stick fluorocarbon O-ring, 5190-2269




Non-stick BTO septa, 5190-3157

Ideal for high usage laboratories, our bulk supplies provide the quality and consistency of Agilent chromatography supplies in convenient and economical packaging. We currently offer Agilent inlet liners, septa, gold inlet seals, and liner O-rings in bulk packaging.

- Economical and convenient packaging
- Overall cost of ownership reduced
- Same great quality Agilent products

Bulk GC Supplies

Description	Unit	Part No.
Single taper split liner, low pressure drop	100/pk	5190-2275
Single taper splitless liner, no wool	100/pk	5190-2270
Single taper splitless liner, glass wool	100/pk	5190-2271
Double taper splitless liner, no wool	100/pk	5190-2272
Certified gold plated seal kit, includes washer	10/pk	5190-2209
Certified non-stick fluorocarbon O-ring	100/pk	5190-2269
Non-stick fluorocarbon O-ring for Flip Top	100/pk	5190-2268
Non-stick BTO septa, 11 mm	400/pk	5190-3157
Non-stick Advanced Green septa, 11 mm	400/pk	5190-3158



Inlet Septa

Septa are available for a variety of different applications and have different upper temperature limits. Lower temperature septa are usually softer, seal better, and can withstand more punctures (injections) than their high-temperature counterparts. If septa are used above their recommended temperatures, they can leak or decompose, causing sample loss, lower column flow, decreased column life and ghosting. To minimize problems:

- Use within the recommended temperature range
- Change regularly
- Install the retainer nut "finger-tight"
- Use septum purge when available
- Use autoinjectors
- Use sharp syringe needles



Premium Non-Stick Septa

Agilent premium non-stick inlet septa are designed and manufactured to provide a reliable non-contaminating seal. Our tri-fold blister pack ensures that each septum remains clean and ready to use.

- Proprietary plasma treatment prevents sticking and unnecessary inlet cleaning
- Innovative blister packaging keeps each septum clean and ready for use
- Center point guides the needle for easy penetration, less coring and longer life
- Precision molding assures accurate fit in the inlet
- Each batch is tested for bleed on Agilent 6890 GC-FID
- Premium formulations selected for sealing and chromatographic cleanliness
- No need to bake septa before using

Summary of Premium Inlet Septum Characteristics

Septum Type	Bleed	Lifetime	Temperature Limits
Non-Stick BTO (Bleed and Temperature Optimized)	◆◆◆	◆	to 400°C injection port temp
Non-Stick Advanced Green	◆◆	◆◆	to 350°C
Non-Stick Long-Life	◆	◆◆◆	to 350°C

◆◆◆ = best ◆◆ = very good ◆ = good



BTO septa, 5183-4757

Non-Stick Bleed and Temperature Optimized (BTO) Septa

- Extended temperature range, lowest bleed
- Maximum injection port temperature 400°C
- Plasma treatment eliminates sticking in the injection port
- Pre-conditioned; ready to use
- Blister packaging for cleanliness and convenience
- Ideal for use with low-bleed, "Mass Spec" capillary columns

Non-Stick Bleed and Temperature Optimized (BTO) Septa

Description	Unit	Part No.
11 mm septa	50/pk	5183-4757
11 mm septa	100/pk	5183-4757-100
11 mm septa	400/pk	5190-3157
5 mm septa through-hole for on-column, in glass jar	50/pk	5183-4758

Comparison of Coring, With and Without CenterGuide (30x magnification)



High-Temperature Septa without CenterGuide

Major coring before 100 autoinjections



Agilent BTO Septa with CenterGuide

Very little coring, even after 700 autoinjections



Advanced green septa, 5183-4759

Non-Stick Advanced Green Septa

- True long-life, high temperature green septa
- More injections per septum
- Plasma treatment eliminates sticking in the injection port
- Maximum injection port temperature 350°C
- High-performance alternative to competitors' "green" septa
- Blister packaging for cleanliness and convenience

Non-Stick Advanced Green Septa

Description	Unit	Part No.
11 mm septa	50/pk	5183-4759
11 mm septa	100/pk	5183-4759-100
11 mm septa	400/pk	5190-3158
5 mm septa through-hole for on-column, in glass jar	50/pk	5183-4760



Long-life septa, 5183-4761




Non-Stick Long-Life Septa

- The preferred septa for autosamplers
- Pre-pierced for extended life and reduced coring
- Ideal for overnight runs
- Up to 400 injections per septum
- Plasma treatment eliminates sticking
- Maximum injection port temperature 350°C
- Soft, 45 durometer, easy on autosampler needles
- Blister packaging for cleanliness and convenience

Non-Stick Long-Life Septa

Description	Unit	Part No.
11 mm septa	50/pk	5183-4761
11 mm septa	100/pk	5183-4761-100
5 mm septa through-hole for on-column, in glass jar	50/pk	5183-4762

Septa Troubleshooting

Symptom	Possible Causes	Remedy
<p>Extra Peaks/Humps</p> 	<p>Septum bleed</p>	<p>Turn off injector heater. If extra peaks disappear, use septum specified for higher temperature or analyze at lower inlet temperature.</p>
<p>Baseline Change After Large Peak</p> 	<p>Large leak at septum during injection and for a short time thereafter (common with large diameter needles)</p>	<p>Replace septum and use smaller diameter needles.</p>
<p>Retention Times Prolonged</p> 	<p>Carrier gas leaks at septum or column connection</p>	<p>Check for leaks. Replace septum or tighten connections if necessary.</p>



General Purpose Septa


Agilent's general purpose septa are made from an enhanced injection-molded silicone rubber. The septa material, gray in color, is specified to withstand over 200 automatic injections at an injection port temperature of 350°C.

General Purpose Septa

Description	Unit	Part No.
Gray Septa		
11 mm septa	50/pk	5080-8896-50
11 mm septa	100/pk	5080-8894-100
9.5 mm (3/8 in.) septa*	50/pk	5080-8728-50
9.5 mm (3/8 in.) septa*	100/pk	5080-8728-100
5 mm through hole septa for on-column inlets, automatic or manual injections**	25/pk	5181-1260
5 mm septa for high column backpressure, on-column inlets**	25/pk	5181-1261

*for 5700 series and 5830/40 GCs

**5 mm septa are packaged in glass jars



Merlin Microseal

- Low bleed, longer life alternative to standard septa for split/splitless injection
- Has a lifetime of more than 2000 injections, depending on samples and operating conditions
- Greatly reduced instrument downtime for septa changes and injection port liner changes due to septa particulates
- Two distinct sealing mechanisms: double O-ring type seal around the syringe needle and spring assisted duckbill to seal the injection port

Merlin Microseal

Description	Part No.
High Pressure Merlin Microseal	
High pressure Merlin Microseal starter kit Includes microseal septum and nut	5182-3442
Merlin Microseal high pressure septum	5182-3444
Microseal high pressure nut	5182-3445
High sample volume septum kit Contains Merlin high pressure Microseal, six 23-gauge syringes, 500 vials and caps	5181-8839
Merlin Microseal Standard Pressure	
Merlin Microseal kit, original low pressure system Includes nut and septum	5181-8816
Merlin Microseal kit, original low pressure system Includes nut and 2 septa	5181-8833
Merlin Microseal septum, stainless steel, rubber (30 psi)	5181-8815
Microseal PTFE nut liners, 2/pk	5182-0853
Merlin Microseal manual syringe, 5 μ L, 23 gauge	5182-3438
Merlin Microseal manual syringe, 10 μ L, 23 gauge	5182-3439
Syringe, 5 μ L, 23 gauge	9301-0892
Syringe, standard plunger, 10 μ L, 23 gauge	9301-0713



For in-depth information about maintaining your GC/MS, request "Maintaining Your Agilent GC and GC/MS Systems" from your Agilent Representative (**publication number 5990-5451EN**).

Inlet Liners

NEW!

Agilent Ultra Inert Liners



Agilent ultra inert liners are the perfect companion to Agilent J&W Ultra Inert GC columns. They provide reproducible inertness liner after liner, maintained through a sequence of samples, and for a range of analytes. Agilent's ultra inert liners were developed – and are manufactured and certified – using a suite of tests specifically designed to ensure batch-to-batch uniformity.

- Exceptional batch-to-batch liner uniformity
- Low to no bleed or background contamination
- Superior coverage, allowing use of glass wool even with highly active compounds

Only ultra inert liners are delivered in Agilent's exclusive touchless packaging with a pre-cleaned, conditioned and non-stick plasma treated O-ring pre-installed. Touchless packaging aids in removal of the old liner, and easy installation of the new, clean, preconditioned liner – without risk of contamination from touching.

Agilent Ultra Inert Liners

Description	Volume (µL)	ID (mm)	1/pk	5/pk	25/pk
Split Inlet Liners					
Straight, ultra inert liner with glass wool	990	4	5190-2294	5190-3164	5190-3168
Splitless Inlet Liners					
Single taper, ultra inert liner	900	4	5190-2292	5190-3162	5190-3166
Single taper, ultra inert liner with glass wool	900	4	5190-2293	5190-3163	5190-3167
Universal Inlet Liners					
Low pressure drop, ultra inert liner with glass wool	870	4	5190-2295	5190-3165	5190-3169

Agilent MS Certified Liners

Agilent MS certified split and splitless liners are manufactured and tested to our highest level of scrutiny to assure reproducibility.

We have built years of experience into MS certified liners to provide the quality and consistency needed for critical applications, especially those using esterification agents for trace level analysis, such as toxicology or drugs of abuse applications.

- Geometrical dimensions and tolerances of the glass are controlled by Statistical Process Control (SPC) with 100% Go-No-Go check
- Glass wool is pre-qualified with mass spectrometry, then inserted using a unique manufacturing procedure to improve reproducibility
- Deactivated MS certified liners are treated with Agilent's proven proprietary deactivation process developed to last longer than other commercially available treatment
- Random samples of MS certified liners are tested using both FID and MSD analysis of challenging probes to evaluate acid/base deactivation, response linearity, peak symmetry, bleed and background noise
- Each Agilent MS certified liner is traceable by the lot codes silk screened on the liner



MS certified single taper split liner, 5188-6576



MS certified straight split liner, 5188-6574



Single taper, glass wool splitless liner,
5062-3587/5188-6568

Agilent MS Certified Liners

Description	Volume (μL)	ID (mm)	1/pk	5/pk	25/pk
Split Inlet Liners					
Single taper, MS certified liner with restriction to hold glass wool	870	4	5188-6576		
Straight, MS certified liner with glass wool	990	4	5188-6574	5188-6569	
Splitless Inlet Liners					
Single taper, MS certified liner with glass wool	900	4	5188-6568	5188-6567	5188-6566

Agilent Split Liners

Agilent single taper split liners are made to strict dimension specifications for optimal inlet performance and feature the tightest tolerances for OD, ID, taper, and glass wool placement. For ease-of-use and reproducibility, the liners have a positioning bead, a restriction to secure the position of the glass wool, and a feature to consistently self-position to the recommended height. The liners also feature Agilent's proprietary deactivation.

Agilent Split Liners




Single taper split liner, 5183-4647, 5183-4711



Straight split liner, 19251-60540



Straight split liner with cup, glass wool, and packing, 18740-60840



Focus liners, 210-4004-5, 210-4022-5

Tips & Tools

Agilent recommends part number 5183-4647 as the top split liner choice for:


- Highest run-to-run area reproducibility
- Least discrimination for wide boiling point range samples
- Use with widest range of conditions and sample types
- Easy self-adjusting installation

Description	Volume (μL)	ID (mm)	1/pk	5/pk	25/pk	100/pk
Single Taper Split Liners						
Single taper, glass wool, deactivated, low pressure drop	870	4	5183-4647	5183-4701	5183-4702	5190-2275
Single taper, MS certified liner with restriction to hold glass wool	870	4	5188-6576			
Single taper, glass wool, deactivated	870	4	5183-4711	5183-4712	5183-4713	
Straight Split Liners						
Straight, glass wool, non-deactivated	990	4	19251-60540	5183-4691	5183-4692	
Straight, MS certified liner with glass wool	990	4	5188-6574	5188-6569		
Straight, with cup (for manual injections)	800	4	18740-80190	5183-4699	5183-4700	
Straight split liner with cup, glass wool, and packing, 18740-60840 (for manual injections, not recommended for use with EPC)	800	4	18740-60840	5183-4697	5183-4698	
Focus Liners						
Deactivated with glass wool	935	4		210-4004-5		
Tapered, deactivated with glass wool	880	4		210-4022-5		
SPME Liners						
Deactivated, for SPME	70	0.75		5188-6471		

Agilent Splitless Liners

Agilent's proprietary deactivation is important for splitless liners because of the longer sample/liner contact time in splitless mode.

Agilent Splitless Liners



Single taper splitless liner, 5181-3316/5181-3316i



Single taper, glass wool splitless liner,
5062-3587/5188-6568



Double taper splitless liner, 5181-3315



Straight, non-deactivated, quartz splitless liner,
18740-80220/5181-8818



Straight, non-deactivated splitless liner, 210-3003



Direct inject liner, 18740-80200


Description	Volume (μ L)	ID (mm)	1/pk	5/pk	25/pk	100/pk
Single Taper Splitless Liners						
Single taper, deactivated	900	4	5181-3316	5183-4695	5183-4696	5190-2270
Single taper, inert	900	4	5181-3316i			
Single taper, glass wool, deactivated	900	4	5062-3587	5183-4693	5183-4694	5190-2271
Single taper, MS certified liner with glass wool	900	4	5188-6568	5188-6567	5188-6566	
Double Taper Splitless Liners						
Double taper, deactivated	800	4	5181-3315	5183-4705	5183-4706	5190-2272
Straight Splitless Liners						
Straight, deactivated, quartz	250	2	5181-8818	5183-4703	5183-4704	
Straight, non-deactivated, quartz	250	2	18740-80220	5183-4707	5183-4708	
Straight, non-deactivated	990	4	210-3003	210-3003-5		
Direct Inlet Liners						
Straight, non-deactivated (for gas samples, headspace, or purge & trap)	140	1.5	18740-80200	5183-4709	5183-4710	

NEW!

Agilent Dimpled Liners

- Single tapered liner
- Agilent proprietary deactivation
- Blocked line of sight to column
- No glass wool
- Optimized for cold splitless or solvent vent
- Recommended for heavy matrix samples, such as pesticides in food extracts

The 2 mm dimpled liner was designed for use in the Agilent MultiMode Inlet for the analysis of heavy matrix samples, such as food safety extracts. Cold splitless injection mode is recommended. Sample residue collects and spreads within the body of the liner – not reaching the taper, bottom of the inlet, or the GC column. This extends the number of analyses before the liner needs to be replaced, improving productivity while reducing inlet maintenance and cost of ownership.



Single taper dimpled splitless liner, 5190-2296

Agilent Dimpled Liners


Description	Volume (µL)	ID (mm)	Part No.
Dimpled splitless single taper, deactivated	200	2	5190-2296

Agilent Direct Connect Liners

Agilent's Direct Connect Liners provide maximum recovery and minimal decomposition of active compounds for methods requiring splitless injection, such as EPA Method 8270. They are best for relatively clean samples containing active analytes, such as water extracts. The liners directly connect with the column, similar to press-fit connectors, to aid complete transfer of sample onto the column, eliminating the problem of inlet discrimination and further increasing sensitivity.



Single taper direct connect liner, G1544-80730



Dual taper direct connect liner, G1544-80700

Direct Connect Liners

Description	Part No.
Single taper direct connect liner, splitless, 4 mm ID, Agilent proprietary deactivation	G1544-80730
Dual taper direct connect liner, splitless, 4 mm ID, Agilent proprietary deactivation	G1544-80700
Single taper direct connect liner, splitless, 4 mm ID, deactivated, inert	G1544-80731



Liner O-rings in dial packaging

Liner O-Rings

- Liners are sealed in the inlet with O-rings or graphite seals
- Graphite seals are used when inlet temperatures exceed 350°C
- Fluorocarbon O-ring seals are easier to replace than graphite that deforms and flakes apart

Only Agilent Fluorocarbon Liner O-rings are:

- Pre-cleaned, then conditioned to eliminate out-gassing of contaminants, which is especially important for trace, ECD and MSD analyses
- Plasma treated for a non-stick, contaminant-free surface that won't stick to the inlet metal surface
- Packaged for convenience and cleanliness in a novel dial package that delivers 1 clean O-ring at a time

Liner O-Rings

Description	Unit	Part No.
Certified non-stick fluorocarbon O-ring	10/pk	5188-5365
	100/pk	5190-2269
Graphite O-ring for splitless liner	10/pk	5180-4173
Graphite O-ring for split liner	10/pk	5180-4168
Non-stick fluorocarbon liner O-ring for Flip Top	10/pk	5188-5366
	100/pk	5190-2268
High temperature PTV inlet liner fluorocarbon O-ring	10/pk	5188-5311



Vespe/Graphite ferrules, 5181-3323

Capillary Column Ferrules and Nuts


Using the wrong ferrule or a worn-out ferrule to seal your column connection can result in inconsistent and unreliable chromatography. An improper ferrule can cause leaks, which allow air and other contaminants to enter the instrument through the column seal, causing major interference with column and detector performance.

For optimum performance, ferrules should be replaced every time the column is replaced and when performing column maintenance.


To minimize problems, follow these general techniques for ferrule installation:

- Don't overtighten — finger tighten the column nut, then use wrench to tighten
- Maintain cleanliness
- Bake out ferrules prior to use (Vespe and Vespe/Graphite only)
- Avoid contamination, such as fingerprint oils
- Inspect used ferrules with magnifier for cracks, chips, or other damage before reusing them
- Change ferrules when new columns or injector/detector parts are installed

Caution: Agilent capillary column ferrules are manufactured to tolerances specific for Agilent J&W HP and DB brand columns. Fused silica tubing used for Agilent J&W VF columns has slightly different tolerances for the outer diameter.



For more information about capillary column compatibility, please visit www.agilent.com/chem/gcsupplies



For the most reliable seals, please refer to the table below for the recommended ferrules for your Agilent GC column.

Ferrule Compatibility

Column ID	DB and HP		VF	
	Graphite	Graphite/Vespe	Graphite	Graphite/Vespe
0.25	500-2114	5181-3323	CR211104	CR213104
0.32	5080-8853	5062-3515	CR211105	CR213105
0.53	5080-8773	5062-3512	CR211108	CR213108


Ferrule Selection Recommendations

Ferrule/Seal Type	Upper Temp. Limit	Usages	Advantages	Limitations
Graphite (100%)	450°C	<ul style="list-style-type: none"> General purpose for capillary columns Suitable for FID and NPD Recommended for high temperature and cool on-column applications 	<ul style="list-style-type: none"> Easy-to-use stable seal Higher temperature limit Can be removed easily 	<ul style="list-style-type: none"> Not for MS or oxygen-sensitive detectors Soft, easily deformed or destroyed Possible system contamination
Vespel/Graphite (85%/15%)	350°C	<ul style="list-style-type: none"> General purpose for capillary columns Recommended for MS and oxygen-sensitive detectors Most reliable leak-free connection 	<ul style="list-style-type: none"> Mechanically robust Long lifetime 	<ul style="list-style-type: none"> Not reusable Flows at elevated temperature Must re-tighten frequently
Vespel (100%)	280°C	<ul style="list-style-type: none"> Isothermal operation Can be reused or removed easily Excellent sealing material when making metal or glass connections 	<ul style="list-style-type: none"> Mechanically robust Long lifetime Can be reused or removed easily 	<ul style="list-style-type: none"> Leaks after temperature cycle Flows at elevated temperature Must re-tighten frequently

Tips & Tools

100% Vespel ferrules should only be used for isothermal applications.

SilTite ferrules are required for leak-tight seals with the Ultimate Union, Deans Switch, and Effluent Splitter.





Graphite ferrules, 5080-8853



Vespel ferrule, 5181-3322



Vespel/Graphite ferrule, 5062-3514

Capillary Column Ferrules and Nuts

Ferrule ID (mm)	Column ID (mm)	Unit	Part No.
General Purpose Graphite Ferrules (Short)			
0.5	0.1, 0.2, 0.25, 0.32	10/pk	5080-8853
0.4	0.05, 0.25	10/pk	500-2114
0.8	0.45, 0.53	10/pk	500-2118
1.0	0.53	10/pk	5080-8773
85% Vespel, 15% Graphite Ferrules (Short)			
0.4	0.1, 0.2, 0.25	10/pk	5181-3323
0.5	0.32	10/pk	5062-3514
0.8	0.45, 0.53	10/pk	5062-3512
Preconditioned 85% Vespel, 15% Graphite Ferrules (Long)*			
0.3	0.1	10/pk	5062-3507
0.4	0.1, 0.2, 0.25	10/pk	5062-3508
0.5	0.32	10/pk	5062-3506
0.8	0.53	10/pk	5062-3538
100% Vespel, High Performance Ferrules (Short)**			
0.4	0.1, 0.2, 0.25	10/pk	5181-3322
0.5	0.32	10/pk	5062-3513
0.8	0.45, 0.53	10/pk	5062-3511
Specialty Ferrules, 85% Vespel, 15% Graphite			
Two Hole			
0.3		10/pk	5181-3388
0.4	0.1, 0.2, 0.25	10/pk	5062-3580
0.5	0.32	10/pk	5062-3581
No hole			
		10/pk	5181-3308
High Temperature PTV Inlet, SS/Graphite			
0.4	0.32	10/pk	5188-5315
0.4	0.53	10/pk	5188-5314


*These ferrules are recommended for use with Agilent GC/MS transfer lines with P/N 05988-20066 MS interface column nut.

**These ferrules are recommended for use in isothermal analysis only.

Tips & Tools

Look for the following signals that indicate ferrule damage:

- Background noise from oxygen diffusing into the system
- Column bleed catalyzed by oxygen
- Sample degradation
- Sample loss
- Increase in detector signal/noise
- Poor retention time reproducibility





SilTite metal ferrules, 5188-5361



Universal column nut, 5181-8830



MS interface column nut, 05988-20066

SilTite Metal Ferrules

Description	Unit	Part No.
For use with 0.25 mm ID capillary columns	10/pk	5188-5361
For use with 0.32 mm ID capillary columns	10/pk	5188-5362
For use with 0.53 mm ID capillary columns	10/pk	5188-5363
For use with 1/16 in. OD stainless steel tubing Includes 2 column nuts	10/pk	5184-3571

Column Nuts

Description	Part No.
Short Nuts	
Universal column nut, 1/16 in. hex, 2/pk	5181-8830
Finger tight column nut for 530 µm columns*	5020-8293
Finger tight column nut for 320 µm columns and smaller*	5020-8292
Blanking plug, finger tight style	5020-8294
6850 column nut, 2/pk	5183-4732
Extended column nut, VI inlet	G3504-20504
High Temperature SimDis PTV inlet, 4 mm hex	5188-5312
Long Nuts	
MS interface column nut, female	05988-20066
Inlet column nut for long or long two-hole ferrules	05921-21170
Accessories	
Open end wrench, 1/4 and 5/16 in.	8710-0510

*For use with graphite ferrules only

Straight Ferrules

Description	Unit	Part No.
1/4 in. PTFE	10/pk	0100-1378
1/4 in. Graphite	10/pk	0100-1324
1/8 in. Graphite	10/pk	0100-1325
1/8 in. 85% Vespel/15% Graphite	10/pk	0100-1332
1/16 in. PTFE	10/pk	0100-1375
1/16 in. Graphite	10/pk	0100-1326
1/16 in. VG-2 Vespel, 40% Graphite	10/pk	0100-1379
6.4 mm Vespel		0100-1104
1/4 in. 85% Vespel/15% Graphite	10/pk	0100-1331

Reducing Ferrules

Description	Unit	Part No.
1/8 to 1/16 in. Vespel	10/pk	0100-1342
1/8 to 1/16 in. VG-1 Vespel, 15% Graphite	10/pk	0100-1344
1/16 in. to 0.4 mm VG-2 Vespel, 40% Graphite	10/pk	0100-1381

Ferrules for LTM Rapid Heating/Cooling System

Description	Unit	Part No.
For use with 0.25-0.4 mm ID LTM columns	5/pk	5190-1437
For use with 0.4-0.5 mm ID LTM columns	5/pk	5190-1438
For use with 0.5-0.8 mm ID LTM columns	5/pk	5190-1439

Ferrules and Nuts for NCD and SCD

Description	Part No.
Spare column nut and ferrule kit	G6600-80018





Capillary Flow Technology Supplies

Agilent offers a family of GC accessories based on our proprietary Capillary Flow Technology. These accessories increase system productivity and performance:

- QuickSwap MS Interface provides vent-free removal of columns
- Deans Switch device simplifies the analysis of complex samples
- Purged Effluent Splitter for inert, leak-free column effluent splitting



Ultimate Union

Ultimate Union

The Ultimate Union is part of Agilent's Capillary Flow Technology family, providing extremely low dead volume column connections. Like the QuickSwap, Deans Switch and Purged Effluent Splitter, the Ultimate Union uses special fittings and SilTite ferrules to create an inert, leak-free and robust seal that doesn't need re-tightening after temperature cycles.

Each Agilent Ultimate Union kit contains:

- 1 Union (your choice of deactivated or non-deactivated)
- 2 Internal nuts
- 1 Swaging nut
- 1 Oven wall clip

Kits do not include SilTite ferrules. Please order ferrules for your column ID separately. SilTite ferrules include 2 transfer line nuts.



Internal nut, G2855-20530



Tee, inert, G3184-60065


Ultimate Union Kits, Fittings and Ferrules

Description	Part No.
Ultimate union kit, deactivated	G3182-61580
Ultimate union kit, non-deactivated	G3182-61581
Internal nut	G2855-20530
Swaging nut	G2855-20555
SilTite metal ferrules, 0.10-0.25 mm ID capillary columns	5188-5361
SilTite metal ferrules, 0.32 mm ID capillary columns	5188-5362
SilTite metal ferrules, 0.53 mm ID capillary columns	5188-5363
Tee, inert	G3184-60065

Tips & Tools

Agilent's QuickSwap Interface Restrictors can increase the productivity of your Agilent 5975 inert MSD system.

Turn to page 326.



Fittings, Ferrules and Supplies

For leak-free, low dead volume and inert column connections with capillary flow accessories, such as the Deans Switch or QuickSwap MS Interface, use only SilTite ferrules and specified nuts. For Capillary Flow devices, use deactivated fused silica tubing. Do not use tubing that has been coated with stationary phase.

Fittings, Ferrules and Supplies

Description	Unit	Part No.
Internal nut		G2855-20530
Swaging nut		G2855-20555
Tee, inert		G3184-60065
Column storage fitting		G2855-20590
SilTite metal ferrules, 0.10-0.25 mm ID capillary columns	10/pk	5188-5361
SilTite metal ferrules, 0.32 mm ID capillary columns	10/pk	5188-5362
SilTite metal ferrules, 0.53 mm ID capillary columns	10/pk	5188-5363
Ferrule pre-swaging tool		G2855-60200

Column/Retention Gap Installation Supplies

Description	Part No.
250 µm retention gap, one 5 m piece	160-2255-5
320 µm retention gap, one 5 m piece	160-2325-5
530 µm retention gap, one 5 m piece	160-2535-5
Fused silica, deactivated, 0.15 mm x 1 m	160-2625-1
Fused silica, deactivated, 0.15 mm x 5 m	160-2625-5
Fused silica, deactivated, 0.15 mm x 10 m	160-2625-10



Glass press-fit connections



Quartz splitter

Press-fit Capillary Column Connectors

In the past it was necessary to use press-fit connectors with specific dimensions to connect columns of those dimensions. Modern press-fit connectors are "laser-milled" to provide highly reproducible taper angles throughout the length of the press-fit, ensuring an excellent seal. Now the only choice you have to make is between a glass union for standard applications, and fused silica unions or deactivated quartz unions for applications demanding maximum inertness.

Glass and Fused Silica Press-fit Connectors

Description	Unit	Part No.
Glass union, universal, 2-way	25/pk	705-0825
Fused silica union, universal, 2-way	5/pk	705-0905
Fused silica union, universal, 2-way	25/pk	705-0925
Fused silica union, universal, 3-way		705-0903
Polymide sealing resin, 5 g		500-1200

Quartz Press-fit Connectors/Splitters

Description	Unit	Part No.
Quartz column connector, 0.1 to 0.53 mm	5/pk	5181-3395
Deactivated quartz column connector	5/pk	5181-3396
Quartz splitter		5181-3397
Quartz deactivated splitter		5181-3398

Mechanical Capillary Column Connectors

Description	Unit	Part No.
Connector, body and nut		5061-5801
Ferrules for Connectors		
Vespel, 0.2 to 0.53 mm ID tubing	2/pk	5061-5804
Vespel, 0.32 to 0.32 mm ID tubing	2/pk	5061-5805
Vespel, 0.32 to 0.53 mm ID tubing	2/pk	5061-5806
Fused silica, undeactivated 530 µm tubing, 10 m		160-2530-10



Graphpak connector for Agilent capillary detectors



Graphpak divider for simultaneous sampling



Capillary injection port connector, 5021-7170

Graphpak Capillary Column Connectors (2.5 mm)*

Column ID (mm)	Connector ID (mm)	Part No.
Capillary Detector Port Connector		
0.32/0.25	0.4	5021-7166
0.53	0.7	5021-7164
Capillary Divider for Simultaneous Sampling		
0.32/0.25	0.53	5021-7148
0.53	0.7	5021-7146
Capillary Injection Port Connector		
0.2	0.3	5021-7169
0.32/0.25	0.4	5021-7170
0.53	0.7	5021-7168

*The 2.5 mm Graphpak is not compatible with the Graphpak 2M used for the PTV.

Note: Order ferrules in addition to the connector to fit your column. Ferrules must be ordered separately.

Ferrules for Connectors

Column ID (mm)	ID (mm)	Unit	Part No.
0.2	0.3	10/pk	5021-7136
0.32/0.25	0.4	10/pk	5021-7137
0.53	0.7	10/pk	5021-7134
Graphpak plug ferrule		10/pk	5021-7133
Replacement Graphpak column nut		5/pk	5062-3525

Valves and Loops



General purpose gas sampling valves

Gas Sampling General Purpose Valves

Description	Part No.
6-port replacement valve WE series, 400 psi, 225°C	5062-9508
6-port replacement valve WE series, Hastelloy C, 400 psi, 225°C	5062-9509
10-port replacement valve WE series, 400 psi, 225°C	5062-9510
10-port replacement valve WE series, Hastelloy C, 400 psi, 225°C	5062-9511
6-port replacement valve WT series, 300 psi, 350°C	0101-0584
10-port replacement valve WT series, 300 psi, 350°C	0101-0585



General purpose liquid sampling valves

Liquid Sampling General Purpose Valves

Description	Part No.
0.2 µL replacement valve UWP series, 1,000 psi, 175°C	0101-0636
0.5 µL replacement valve UWP series, 1,000 psi, 175°C	0101-0637
1.0 µL replacement valve UWP series, 1,000 psi, 175°C	0101-0638
0.5 µL replacement valve UW series, 5,000 psi, 75°C	0101-0639




Front ferrules, stainless steel, 5181-1292

Parts for Interfacing Capillary Columns to W-Series Valves

Description	Part No.
Stainless steel bulkhead ZDV union, for interfacing 530 μ m columns with 1/16 in. transfer line	0100-1515
Stainless steel bulkhead ZDV union, for interfacing 320 μ m columns with 1/16 in. transfer line	0100-1527
Polymide ferrule, 1/16 in.	0100-1512
Polyimide liner for 530 μ m columns	0100-1513
Polyimide liner for 320 μ m columns	0100-1514
1/16 in. stainless steel counterbored nut	0100-1511
Installation tool for liners	18900-20850

Valve Supplies

Description	Part No.
1/16 in. stainless steel nut	5181-1291
1/16 in. front ferrule, stainless steel	5181-1292
Straight metering valve, 1/16 in., stainless steel, for LSVs as a sample-out restrictor or as a flow-balancer for 10-100 mL/min	0101-0355
Micrometering valve, for flow balancing gas flows of 2-50 mL/min	0101-0633
Air actuator	19325-60660
Solenoid valve for controlling actuator	05890-61090
Solenoid valve for controlling actuator, for GCs after serial number 3223A43573	05890-61095
Angle metering valve, 1/16 in., stainless steel	0101-0403




Valve Loops for GC (Includes loop, nut and ferrule, 1/16 in.)

Description	Part No.
Sample loop, 0.25 cc	0101-0303
Sample loop, 0.50 cc	0101-0282
Sample loop, 1.00 cc	0101-0299
Sample loop, 2.00 cc	0101-0300
Sample loop, 5.00 cc	0101-0301
Sample loop, 10.00 cc	0101-0302
A. Tube, 1/16 in. stainless steel, 280 mm long	18900-20250
Tube, 1/16 in. stainless steel, 400 mm long	18900-20280
Tube, 1/16 in. stainless steel, 375 mm long	18900-20281
Tube, 1/16 in. stainless steel, 560 mm long	18900-20300
B. Tube, 1/16 in. with 1/8 in. flare, stainless steel, 360 mm long	1530-2163
Tube, 1/16 in. with 1/8 in. flare, stainless steel, 520 mm long	1530-2167
C. Tube, 1/16 in. with 1/8 in. bulkhead fitting, 520 mm long	07675-80050
Tube, 1/16 in. with 1/8 in. bulkhead fitting, nickel, 460 mm long	18900-80255

Replacement Rotors for Gas Sampling Valves

Description	Part No.
6-port replacement rotor WE series, 400 psi, 225°C	5181-7459
10-port replacement rotor WE series, 400 psi, 225°C	5181-7460
6-port replacement rotor WT series, 300 psi, 350°C	1535-4952
10-port replacement rotor WT series, 300 psi, 350°C	1535-4954

Sample Introduction Systems



7693A Automatic Liquid Sampler

7693A Automatic Liquid Sampler Replacement Parts and Supplies

To support the higher productivity, performance, and flexibility offered by the 7693A ALS, Agilent has expanded its supplies offering. Agilent Blue Line autosampler syringes are specifically designed to support the 7693A, while increasing plunger life and reducing costly downtime. For cost-conscious laboratories, economical shell vials and caps provide quality at an attractive price. Additional accessories, such as color-coded sample trays and vial caps, add to system ease-of-use.


7693A Replacement Parts and Supplies

Description	Unit	Part No.
Gripper finger caps	16/pk	G4514-60710
Injector mounting post		G4513-20561
Dual parking post for autosampler		05890-61525
Needle support insert, standard		G4513-40525
Needle support insert, on-column		G4513-40529
Vial rack, set of 3 Includes 12 label tags (4 colors)		G4514-67505
Vial rack label kit		G4525-60701
Vial rack label kit, red	3/pk	G4525-60702
Vial rack label kit, yellow	3/pk	G4525-60703
Vial rack label kit, green	3/pk	G4525-60704

Tips & Tools

Agilent Blue Line autosampler syringes are designed to support the 7693A.

Turn to page 64.



Shell Vials for 7693A Automatic Liquid Sampler

Description	Unit	Part No.
Clear, glass shell vials and caps	200/pk	5190-1570

Wash Vials (also for standards, diluents)

Description	Unit	Part No.
4 mL wash vials with fill markings and caps	25/pk	5182-0551
Diffusion caps for 4 mL vials	12/pk	07673-40180
Septa for 4 mL vials*	144/pk	9301-1031


*Septa for 4 mL vials should only be used for sample storage



Automatic Liquid Sampler Supplies

Automatic Liquid Sampler Supplies

Description	Unit	Part No.
4 mL wash vial with screw caps	144/pk	9301-0723
Diffusion caps for 4 mL vials	12/pk	07673-40180
Septa for 4 mL vials	144/pk	9301-1031
4 mL wash vials with fill markings and caps	25/pk	5182-0551
Screw for mounting syringe		07673-20570
Quadrant tray (4 tray sections)		18596-40015
7673 Basic Supply Kit		07673-60840
Contains 10 μ L syringes (6/ea), 23/26 gauge needles, 4 mL vials with diffusion caps (144/pk), 2 mL automatic sampler vials with screw caps (1,000/pk), GC septa (25/pk), vial racks (5/pk)		



Bar code reader labels

Bar Code Reader Labels

Description	Part No.
Labels numbered (1,000/roll)	
1 to 1,000	5958-9450
1,001 to 2,000	5958-9441
2,001 to 3,000	5958-9442
3,001 to 4,000	5958-9443
4,001 to 5,000	5958-9444
5,001 to 6,000	5958-9445



7697A Headspace Sampler

7697A Headspace Sampler Supplies


The new 7697A Headspace Sampler from Agilent utilizes advanced designs based on our industry-leading gas chromatography architecture. The headspace sampling technique allows introduction of volatile compounds to the GC or GC/MS from virtually any sample matrix, while leaving unwanted components in a disposable sample vial. With up to 111 sample vial positions and removable vial racks, the 7697A supports nearly continuous operation to satisfy even the busiest laboratory.

- Built-in legendary Agilent pneumatics for superior control and easier setup
- Proven valve and loop sampling technology
- Fully-automatic sample vial leak checking and available bar code reader help ensure greater confidence in results method compatibility
- Instrument control software that is fully integrated in Agilent data systems
- Resource conserving programmable instrument scheduler

7697A Headspace Replacement Parts and Supplies

Description	Part No.
Tray vial racks	G4556-60019
Vial rack label	G4556-90500
Universal/external split vent trap with 3 cartridges, 1/8 in. Swagelok fitting	RDT-1020
Leak test kit	G4556-67010
Includes instruction sheet, no-hole ferrule, 1/8 in. nylon tube fitting plug, headspace leak test vial, 1/16 in. stainless steel ZDV plug, 11 mm low bleed septa (5/pk)	
Sample probe, deactivated SN 1030	G4556-60690
Sample probe, deactivated SN 2000	G4556-60125
6-port valve, replacement rotor, WT series, 300 psi, 350°C	1535-4952

(Continued)





7697A Headspace Sampler

The transfer line heater assembly is 1 m in length and accommodates the following tubing types:

- Fused silica capillary of 0.25 mm, 0.32 mm, and 0.53 mm ID with maximum OD of 0.67 mm
- Metal capillary of 0.53 mm ID, such as Agilent UltiMetal or ProSteel, with maximum OD of 0.67 mm

For one transfer line, a piece of fused silica or ProSteel approximately 1 m in length is required in addition to one ferrule and one nut and reducing union. Order a ProSteel sleeve to protect the transfer line when operating above 200°C. ProSteel operated above 200°C in the transfer line without the sleeve can permanently bind to the heated conduit tube.

7697A Headspace Replacement Parts and Supplies

Description	Part No.
Transfer Line Components	
Deactivated fused silica, 5 m length	
0.25 mm	160-2255-5
0.32 mm	160-2325-5
0.45 mm	160-2455-5
0.53 mm	160-2535-5
ProSteel deactivated stainless steel, 5 m length	
0.53 mm	160-4535-5
Polyimide sleeve for ProSteel	4177-0607
Polyimide, Valcon ferrule, 5/pk	
1/32 in. for tubing OD 0.50 ≤ 0.80 mm	0100-2595
1/32 in. for tubing OD 0.25 ≤ 0.40 mm	0100-2610
Nut and reducing union for 6 port valve and transfer line connection	0100-2594
Septum nut, transfer line, split/splitless and multimode inlets	G3452-60835
Sampling Loops, SN 2000	
Sample loop, 0.025 mL	G4556-80101
Sample loop, 0.05 mL	G4556-80102
Sample loop, 0.10 mL	G4556-80103
Sample loop, 0.50 mL	G4556-80105
Sample loop, 1.00 mL	G4556-80106
Certified sample loop, 1.00 mL	G4556-80126
Sample loop, 3.00 mL	G4556-80108
Certified sample loop, 3.00 mL	G4556-80128
Sample loop, 5.00 mL	G4556-80109
Sampling Loops, SN 1030	
Sample loop, 0.025 mL	G4556-80111
Sample loop, 0.05 mL	G4556-80112
Sample loop, 0.10 mL	G4556-80113
Sample loop, 0.50 mL	G4556-80115
Sample loop, 1.00 mL	G4556-80116
Sample loop, 3.00 mL	G4556-80118
Sample loop, 5.00 mL	G4556-80119

Vials and Caps for 7697A

Description	Unit	Part No.
Vial kit	100/pk	5182-0840
20 mL Headspace crimp top, flat bottom vials, silver aluminum one-piece crimp caps with safety feature, PTFE/white silicone septa		
Headspace crimp top, flat bottom vials, 10 mL	100/pk	5182-0838
Headspace crimp top, flat bottom vials, 20 mL	100/pk	5182-0837
Silver aluminum crimp caps with 20 mm, PTFE/silicone septa	100/pk	5183-4477
20 mm electronic crimper		5062-0208
Ergonomic manual crimper for 20 mm caps		5040-4669




G1888A Headspace unit

G1888A Network Headspace Sampler Supplies

Description	Part No.
Stainless Steel Sample Loops	
Certified sample loop, 1 mL, deactivated	5190-2265
Certified sample loop, 3 mL, deactivated	5190-2266
Sample loop, 1 mL, deactivated	2321700003
Sample loop, 3 mL, deactivated	2321700004
Probes and Unions	
Sample probe, deactivated	2322700011
M6 union, brass	2302533140
Union, zero dead volume, deactivated	2307230001
Union	2307232901
Transfer Line Needles and Unions	
Needle only, headspace transfer line, deactivated 0.5 mm OD	2322590004
Needle only, headspace transfer line, deactivated 0.7 mm OD	2322590005
Strain relief septum nut	6410090050
Tubing	
Tubing, solenoids to 6-port valve, deactivated	0410105017
Tubing, probe to 6-port valve, deactivated	1300502506
Standards	
OQ/PV Headspace Sample	5182-9733
Contains 0.2-0.3% t-butyl disulfide, 1,2-dichlorobenzene, and nitrobenzene in ethanol	
PM Kits	
G1888A PM kit with 1 mL loop	G1888-60702
G1888A PM kit with 3 mL loop	G1888-60703
G1888A enhanced PM kit	G1888-60704

Tips & Tools

For a complete selection of headspace vials, turn to page 49.

**7694 Headspace Sampler Supplies**

Description	Part No.
Needles	
Needle only, headspace transfer line, deactivated 0.5 mm OD	2322590004
Needle for transfer line, 0.25 mm ID, 0.5 mm OD, nickel	301-016-HSP
Needle only, headspace transfer line, deactivated 0.7 mm OD	2322590005
Needle for transfer line, 0.4 mm ID, 0.8 mm OD, nickel	301-015-HSP
Needle assembly vial probe, deactivated	301-220-HSP
Needle assembly (vial probe)	301-013-HSP
Fittings	
Zero dead volume union	325-045-HSP
Transfer line nut	19258-20830
Transfer line ferrule	19258-20870
Union FF 6MB, 5-piece set	325-062-HSP
Union T6 MB, 5-piece set, brass	325-132-HSP
Union T5 MA	325-185-HSP
Valves	
Valve, solenoid vent Kalrez	3600500001
Valve, solenoid vial pressurization	3600500002
Tubing and Transfer Lines	
Sample loop, 1 mL, deactivated	2321700003
Sample loop, 1 mL, nickel	321-055-HSP
Sample loop, 3 mL, deactivated	2321700004
Sample loop, 3 mL, nickel	321-056-HSP
Tube, needle, 6-port valve, deactivated	301-212-HSP
Tube, needle, 6-port valve, nickel	301-169-HSP
Tube, vent-valve stainless steel	301-170-HSP
Sensor tube, 125 mm PTFE	321-057-HSP
Restrictor, stainless steel	321-002-HSP
Transfer line, deactivated, 1 m	301-211-HSP
Transfer line, 1 m, nickel	301-152-HSP
Transfer line, 80 cm, nickel	301-011-HSP
Repair, Leak Test, and OQ/PV Supplies	
Strain relief septum nut	301-205-HSP
Oven adaptor for 10 mL vials	301-017-HSP
Tray adaptors for 10 mL vials, 25/pk	300-301-HSP
Headspace leak test kit	G1888-60701
OQ/PV Headspace Sample	5182-9733

G1883A Network Headspace Supplies

Description	Part No.
Needles	
Needle only, headspace transfer line, deactivated 0.5 mm OD	2322590004
Needle for transfer line, 0.25 mm ID, 0.5 mm OD, nickel	301-016-HSP
Needle only, headspace transfer line, deactivated 0.7 mm OD	2322590005
Needle for transfer line, 0.4 mm ID, 0.8 mm OD, nickel	301-015-HSP
Needle assembly vial probe, deactivated	232-2790012-EHS
Needle assembly vial probe, nickel	232-2790010-EHS
Fittings	
Union elbow M5	998-0000053-EHS
Transfer line nut	19258-20830
Transfer line ferrule	19258-20870
Union FF 6MB, 5-piece set	325-062-HSP
Union T6 MB, 5-piece set, brass	325-132-HSP
Union T5 MA	325-185-HSP
Valves	
Restrictor, stainless steel	321-002-HSP
Valve, solenoid vent Kalrez	3600500001
Valve, solenoid vial pressurization	3600500002
Tubing and Transfer Lines	
Sample loop, 1 mL, deactivated	2321700003
Sample loop, 1 mL, nickel	321-055-HSP
Sample loop, 3 mL, deactivated	2321700004
Sample loop, 3 mL, nickel	321-056-HSP
Oven adaptor for 10 mL vials	301-017-HSP
Tube, needle, 6-port valve, deactivated	301-212-HSP
Tube, needle, 6-port valve, nickel	301-169-HSP
Tube, vent-valve stainless steel	301-170-HSP
Sensor tube, 125 mm PTFE	321-057-HSP
Transfer line, deactivated, 1 m	301-211-HSP
Transfer line, 1 m, nickel	301-152-HSP
Transfer line, 80 cm, nickel	301-011-HSP
Repair, Leak Test, and OQ/PV Supplies	
Strain relief septum nut	301-205-HSP
Headspace leak test kit	G1888-60701
OQ/PV Headspace Sample	5182-9733




Stratum PTC Sample Concentrator



Trap, BTEX + MTBE, 5188-8813

Tips & Tools

Compared to a frit sparger, the fritless sparger may be the better choice when a water sample has a tendency to foam. This sparger is not appropriate for soil samples, which tend to clog the capillary tube.



Teledyne Tekmar Purge and Trap Supplies

Glassware for Teledyne Tekmar Purge and Trap Concentrators, 1/2 in. Mount

Description	Part No.
5 mL frit sparger (glassware only)	5182-0852
5 mL frit sparger kit with fittings	5182-0846
25 mL frit sparger (glassware only)	5182-0851
25 mL frit sparger kit with fittings	5182-0845
5 mL fritless sparger (glassware only)	5182-0850
5 mL fritless sparger kit with fittings	5182-0844
25 mL fritless sparger (glassware only)	5182-0849
25 mL fritless sparger kit with fittings	5182-0796
5 mL needle sparger (glassware only)	5182-0848
5 mL needle sparger kit	5182-0795
25 mL needle sparger (glassware only)	5182-0847
25 mL needle sparger kit	5182-0794

Traps for Teledyne Tekmar Stratum Purge and Trap Concentrator

Description	Part No.
Trap, BTEX + MTBE	5188-8813
Trap #5, OV-1/Tenax/Silica Gel/Charcoal	5188-8814
Trap #8, Carbowax B/Carbosieve S-III	5188-8815
Trap #9, Proprietary	5188-8816
Trap, Tenax/Silica Gel/Carbosieve S-III	5188-8817
Strat-Trap, Tenax/Silica Gel, #2	5188-8818
Strat-Trap, Tenax/Silica Gel/Charcoal, #3	5188-8819
Strat-Trap, OV-1/Tenax, #7	5190-1445
Strat-Trap, Tenax, #1	5190-1446
Trap, VOCARB 3000	5188-8820
Trap, VOCARB 4000	5188-8821
Trap, BTEX	5188-8822

Stratum traps are U-shaped

Atomx VOC Autosampler Supplies

Description	Part No.
Antifoam agent, Antifoam 1520, 10 mL	5190-2235
Syringe with side port, 27 mL	5190-2234
Vessel, amber IS, 15 mL	5190-2233
Frit sparge glassware kit, 25 mL	5190-2232
Fritless sparge glassware kit, 25 mL	5190-2231

Traps for Teledyne Tekmar Velocity Purge and Trap Concentrator

Description	Part No.
Trap, Vocab 3000 (K Trap)	5182-0775
Trap, Vocab 4000 (I Trap)	5182-0774
Trap, Tenax (A Trap)	5182-0783
Trap, Tenax/Silica Gel/Charcoal (C Trap)	5182-0781
Trap, BTEX	5182-0773
DryFlow moisture trap	14-8911-003


Velocity traps are straight



Markes Thermal Desorption system



Sampling tube, MKI-UTD-5105



Silcosteel difflok cap, MKI-MTD-1204

Markes Thermal Desorption

Agilent now offers a comprehensive line of supplies for Markes Thermal Desorption (TD) instrumentation. Thermal desorption allows the introduction of volatile and semi-volatile compounds from a wide range of sample matrices, directly into a GC or GC/MS.

Markes Thermal Desorption Supplies

Description	Unit	Part No.
O-rings, Markes 7 mm cold trap seals	10/pk	MKI-U-COV07
O-rings, Markes 6 mm cold trap seals	10/pk	MKI-U-COV06
PTFE filter disks, 5.1 mm Markes TD	10/pk	MKI-U-DISK1
PTFE filter disks, 6.3 mm Markes TD	10/pk	MKI-U-DISK3
Spare general purpose carbon cold trap		MKI-U-T11GPC
Sampling tube, Tenax TA, Markes Unity		MKI-UTD-5105
Quick fit connectors, Markes Unity	10/pk	MKI-C-QSC10
Stainless steel Difflok cap, Markes Unity	10/pk	MKI-MTD-1169
Silcosteel Difflok cap, Markes Unity	10/pk	MKI-MTD-1204
O-ring insertion tool, Markes Unity TDI		MKI-Z-0285
O-ring extraction tool, Markes Unity TDI		MKI-Z-0351
Cold trap alignment tool, Markes Unity		MKI-UTD-5064
Cold trap, air toxics, C ₂ -C ₁₄ , Unity 2		MKI-U-T3ATX-2S
Cold trap, air toxics, C ₂ -C ₁₄ , Unity		MKI-U-T3ATX
Cold trap, materials emissions, Unity		MKI-U-T12ME
Cold trap, GP Carbon, C _{4/5} -C _{30/32} , Unity 2		MKI-U-T11GPC-2S
O-rings, 010 Markes Unity	10/pk	MKI-U-COV10
Cold trap, materials emissions, Unity 2		MKI-U-T12ME-2S
Empty stainless steel TD tubes	10/pk	C-TBE10
Tenax stainless steel tubes, preconditioned/capped	10/pk	C-TBP1TC
Empty glass TD tubes	10/pk	C-GT010
PTFE inserts	10/pk	C-PL010
Long term TD tube storage caps	10/pk	C-CF020
Cap-LOK Tool for long term storage caps		C-CPLOC
Diffusive sampling caps	10/pk	C-DF010
Bio-VOC breath samplers	10/pk	C-BIO10
Disposable card mouth piece for Bio-VOC	10/pk	C-B010M
Tenax TA 34-60 Mesh, 10 g		C-TNXTA

(Continued)



Markes Thermal Desorption Supplies

Description	Unit	Part No.
General purpose hydrophobic tubes, stainless steel Preconditioned and capped with 1/4 in. brass storage caps. For pumped sampling n-C ₅ to n-C ₂₀ .	10/pk	C-HY010C
Tenax/S'carb 'Sulphur' tubes Preconditioned and capped with 1/4 in. brass storage caps. For odor and landfill gas analysis.	10/pk	C-102SSC
Carbograph 1 stainless steel tubes Preconditioned and capped with 1/4 in. brass storage caps. For pumped sampling C ₅ to C ₁₄ , plus diffusion of BTX.	10/pk	C-TBP1C1C
Carb X stainless steel tubes Preconditioned and capped with 1/4 in. brass storage caps. For pumped/diffusion of 1,3-butadiene & benzene.	10/pk	C-TBP1CXC
Air toxics (TO-17) stainless steel tubes Preconditioned and capped with 1/4 in. brass storage caps. For pumped sampling VOCs n-C ₃ to n-C ₁₂ .	10/pk	C-AT010C
Universal stainless steel tubes Preconditioned and capped with 1/4 in. brass storage caps	10/pk	C-UN010C
Glass tubes with 1 cm Tenax For direct liquid injection	10/pk	C-G1CM10
Glass air toxics (TO-17) tubes Pre-packed with 2 carbon-based sorbents; preconditioned and capped with 1/4 in. brass storage caps	10/pk	C-GAT010C
CRS BTX Standards, 1 µg	10/pk	C-BTX1UG

Inlet Systems

7890 Turn Top Inlet System




Convenient new turn top design is built into every 7890 Split/Splitless and Multimode Inlet, allowing you to change liners in less than 30 seconds without special tools or training.


7890 Turn Top Inlet System

Description	Unit	Part No.
Turn top		G3430-40035
Split ring		0510-1306
Certified non-stick fluorocarbon O-ring	10/pk	5188-5365
	100/pk	5190-2769

Tips & Tools

To learn more about Agilent's complete portfolio of services, please visit www.agilent.com/chem/services





Flip Top Inlet Sealing System

Agilent's Flip Top Inlet Sealing System is the faster, smarter way to change inlet liners on Agilent 6890, 6850 and 5890 GC systems.

- Cuts liner replacement time to as little as 30 seconds
- Eliminates frustrating searches for special wrenches or tools
- Improves inlet ergonomics – no more handling of heated parts, no more burns or scrapes
- Decreases downtime and increases productivity
- Minimizes exposure to ambient air, extending column life
- Easily installed by user in 15 minutes


Available exclusively from Agilent, the Flip Top has a levered arm that attaches to any 6890/6850/5890 insert weldment and locks to the injection port using an adapter ring screwed onto the inlet. Once installed, simply lift the arm of the Flip Top which releases the insert weldment from the injection port, and allows instant access to the liner. The process is simply reversed to reseal the weldment to the port.



Flip Top Inlet Sealing System installation kit,
5188-2717

Flip Top Inlet Sealing System

Description	Unit	Part No.
Flip Top Inlet Sealing System For 6890, 6850, 5890 only; not compatible with 7890		5188-2717
Non-stick fluorocarbon liner O-ring for Flip Top	10/pk	5188-5366
	100/pk	5190-2268



QuickPick Splitless PM Kit, 5188-6497



QuickPick Purged Packed PM Kit, 5188-6498



Split Vent Trap PM Kit, 5188-6495

Agilent Inlet Convenience Kits

Convenience kits are an easy way to get all the supplies you need using one part number. Agilent's new PM kits include septa, liners, O-rings, gold seals and traps.

Agilent Inlet Convenience Kits

Description	Part No.
QuickPick Split Inlet PM Kit Includes 5 non-stick BTO septa, 1 split liner, 1 non-stick liner O-ring, and inlet gold seal kit	5188-6493
QuickPick Split Vent and Inlet PM Kit Includes 5 non-stick BTO septa, 1 split liner, 1 non-stick liner O-ring, inlet gold seal kit, and split vent trap with 2 O-rings	5188-6496
QuickPick Splitless Inlet PM Kit Includes 5 non-stick BTO septa, 1 splitless liner, 1 non-stick liner O-ring, and inlet gold seal kit	5188-6494
QuickPick Splitless Vent and Inlet PM Kit Includes 5 non-stick BTO septa, 1 splitless liner, 1 non-stick liner O-ring, inlet gold seal kit, and split vent trap with 2 O-rings	5188-6497
QuickPick Purged Packed Inlet PM Kit Includes 5 non-stick BTO septa, 1 O-ring, 1 ferrule, and 1 disposable glass liner	5188-6498
Internal split vent trap PM kit for split/splitless, volatiles, and PTV septumless inlet split vent line Includes 1 cartridge and 2 O-rings	5188-6495

Split/Splitless Inlet Seals

To ensure that you have a consistent and inert surface to properly seal the inlet and prevent sample degradation, Agilent has revolutionized production of the gold inlet seal. Unlike traditional machined seals, the new format Agilent Gold Inlet Seal has a very reproducible smooth surface, eliminating traces of machining grooves that can be the source of minute leaks. With Agilent's proprietary metal injection molding (MIM) manufacturing process, every gold inlet seal provides a high quality, leak-free seal so critical for reproducible results.

Our new package keeps the gold seal clean and scratch-free. For your added convenience, an inlet washer is provided with each inlet seal.



Gold plated seal kit, 5188-5367


Split/Splitless Inlet Seals

Description	Unit	Part No.
Certified gold plated seal kit, includes washer		5188-5367
	10/pk	5190-2209
Gold plated seal with cross*		5182-9652
Stainless steel seal		18740-20880

*Use with total flow rates above 200 mL/min

Split/Splitless Inlets

The combined split/splitless inlet is the most popular inlet for capillary column gas chromatography. Because it can be used in either split or splitless mode, it provides a very effective combination that can cover most analysis requirements.



Split/Splitless Inlet assembly

7890/6890/6850 Split/Splitless Inlet Supplies

Item	Description	Unit	Part No.
	QuickPick Split Inlet PM Kit		5188-6493
	QuickPick Split Vent and Inlet PM Kit		5188-6496
	QuickPick Splitless Inlet PM Kit		5188-6494
	QuickPick Splitless Vent and Inlet PM Kit		5188-6497
1	Headspace septum retainer nut		18740-60830
	Septum retainer nut		18740-60835
	Strain relief septum nut		301-205-HSP
2	11 mm Certified BTO septa	50/pk	5183-4757
	For complete offering of premium septa, see pages [4177]		
3	7890 Insert Weldment		
	Top insert weldment assembly, standard		G3452-60730
	Top insert weldment assembly, headspace		G3452-60100
	Top insert, AC gang fitting weldment		G3430-60011
	Top insert assembly, valve		G3480-67585
6890 Insert Weldment	S/SL insert weldment		G1544-60585
	Used with large charcoal canister type filter, for 6890/6850		
	S/SL insert assembly for G1540A with valved system option		G1580-60585
	This insert assembly uses the large charcoal canister split vent filter, for 6890/6850		
	Similar to G1544-60575 except carrier lines separated for interface to valved systems of a G1540A instrument		G1580-60575
	Original standard EPC using 1/4 in. split vent filter		G1544-60575
	Similar to G1544-60575 except allows insertion of 1/4 in. chemical filters to clean carrier gas for ECD operation		G1544-80580
	Insert Weldment Standard manual pneumatics		19251-60575

(Continued)

7890/6890/6850 Split/Splitless Inlet Supplies

Item	Description	Unit	Part No.
4	Certified non-stick fluorocarbon O-ring	10/pk	5188-5365
	Graphite O-ring for split liner	10/pk	5180-4168
	Graphite O-ring for splitless liner	10/pk	5180-4173
5	Split liner, single taper, low pressure drop, glass wool	1/pk	5183-4647
		25/pk	5183-4702
	Splitless liner, single taper, without glass wool	1/pk	5181-3316
		25/pk	5183-4696
	For complete offering of liners, see pages [4173]		
6	Split vent trap kit		G1544-60610
	Replacement cartridge		G1544-80530
	Includes 2 cartridges and 4 O-rings		
	Split vent trap PM kit		5188-6495
	Includes 1 cartridge and 2 O-rings		
7	Retaining nut		G1544-20590
8	Stainless steel seal		18740-20880
	Certified gold plated seal kit, includes washer ¹		5188-5367
	Replacement for 18740-20885		
	Gold plated seal with cross ²		5182-9652
9	Washers, 0.375 OD	12/pk	5061-5869
10	Reducing nut		18740-20800
11	Insulation kit, 3 pieces		5188-5241
12	Lower insulation cover		19243-00070
13	Ferrules		
	For complete offering of ferrules, see pages [8580]		
14	Universal column nut	2/pk	5181-8830
	6850 column nut	2/pk	5183-4732
	Split/splitless septum nut angled wrench		19251-00100
	Flip Top Inlet Sealing System		5188-2717
	For 6890, 6850, 5890 only; not compatible with 7890		
	Capillary Inlet Evaluation Sample (Split Mode)		8500-4789

¹Use with total inlet flow rates below 200 mL/min²Use with total flow rates above 200 mL/min

(Continued)



Item 6, split vent trap, G1544-60610



Reducing nut, 18740-20800




Gold plated seal kit, 5188-5367

7890/6890/6850 Split/Splitless Inlet Supplies

Item	Description	Unit	Part No.
	Capillary Inlet Supplies Kit, Includes:		5181-8838
	Certified gold plated seal kit, includes washer		5188-5367
	Liner, split, straight, glass wool, non-deactivated	5 each*	19251-60540
	Liner, splitless, single-taper, glass wool, deactivated	2 each*	5062-3587
	Certified non-stick fluorocarbon O-ring	10/pk*	5188-5365
	Liner, direct, 2 mm ID, deactivated		5181-8818
	11 mm Certified BTO septa	50/pk*	5183-4757
	Capillary inlet cleaning wires	5/pk*	5180-4153

*Quantity when part ordered individually



Gold seal on the split/splitless inlet


Multimode Inlet

Agilent's premium inlet – two inlets in one
for maximum performance and flexibility for the 7890A

Whether you need to increase your system sensitivity with large volume injection capabilities, analyze thermally labile compounds or inject dirty samples, the Multimode Inlet delivers the performance and flexibility you need.

Tips & Tools

The multimode inlet and the split/splitless inlet use the same supplies including septa, liners, and O-rings. Please refer to split/splitless inlet ordering information for more details. Turn to pages 296–298.




Exploded Parts View of the Multimode Inlet

Cool On-Column Inlets

Cool on-column injection is superior in many ways to other sample introduction techniques.

Advantages include:

- Elimination of sample discrimination
- Elimination of sample alteration
- Solvent focusing of early eluting solutes
- High analytical precision



Cool On-Column Inlet assembly

7890/6890 Cool On-Column Inlet Supplies

Item	Description	Unit	Part No.
Automatic Injection			
1a	Septum nut for 320 μ m columns		19245-80521
1b	Septum nut base for 530 μ m assembly		G1545-80520
2	5 mm Advanced Green septa	50/pk	5183-4760
	5 mm BTO septa	50/pk	5183-4758
Manual Injection			
1c	Cooling tower assembly		19320-80625
1d	Duck bill	10/pk	19245-40050
	Fused silica syringe needles	6/pk	19091-63000
	On-column syringe, fused silica (barrel only)		9301-0658
Common Supplies			
3	Spring		19245-60760
4	Inserts for capillary columns		
	For 200 μ m columns, 1 ring		19245-20510
	For 250 μ m columns, 6 rings		19245-20515
	For 320 μ m columns, 5 rings		19245-20525
	For 530 μ m columns, no rings		19245-20580
	For 530 μ m Al clad columns, 4 rings		19245-20780
5	320 μ m, 0.5 mm ID graphite ferrule		5080-8853
6	Universal column nut	2/pk	5181-8830


Programmed Temperature Vaporizer (PTV) Inlets

PTV inlets combine the benefits of split, splitless and on-column inlets. The sample is usually injected into a cool liner, so syringe needle discrimination does not occur. Then the inlet temperature is increased to vaporize the sample. The user programs vent times and temperature to achieve the equivalent of split or splitless transfer of sample vapors to the column. PTV injection is considered the most universal sample introduction system because of its flexibility.

7890/6890 Septumless PTV Inlet Supplies

Description	Column ID (mm)	Unit	Part No.
Microseal high pressure nut			5182-3445
Merlin Microseal			5182-3444
Septumless head			G2617-60507
Septum head			G2618-80500
Septum retainer nut			18740-60835
PTV inlet assembly			G2617-60506
PTV LC ₂ cooling jacket			G2617-60508
PTV LN ₂ cooling jacket			G2619-60501
Silver seal		5/pk	5182-9763
Graphpak 2M inlet adapter	0.20		5182-9754
	0.25-0.33		5182-9761
	0.53		5182-9762
Ferrules for Graphpak 2M inlet	0.20		5182-9756
	0.25		5182-9768
	0.32		5182-9769
	0.53		5182-9770
Replacement Graphpak column nut			5062-3525
PTV insulation block			G2617-20510
PTV Cryo insulator			G2617-60510
Teflon ferrule (needle seal)		10/pk	5182-9748
Kalrez seal			5182-9759
Valve body			5182-9757
Pressure spring			5182-9758
Viton seal		5/pk	5182-9775
Sealing element			5182-9760
CO ₂ Cryo inline filter			3150-0602
Service kit for septumless head			5182-9747
Contains Kalrez seal, valve body, and pressure spring			
Graphpak 3D ferrules		5/pk	5182-9749
Assembly tool for Graphpak 3D ferrules			G2617-80540

Septumless Programmable Temperature Vaporizing (PTV) Inlet assembly



Programmable Temperature Vaporizing (PTV) Liners



Description	ID (mm)	Volume (μL)	Part No.
Liners for Septumless PTV Inlet, G3501A, G3502A, G3503A			
PTV liner, single baffle, glass wool, deactivated	2	180	5183-2038
PTV liner, single baffle, deactivated	2	200	5183-2036
PTV liner, multi baffled, deactivated	1.8	150	5183-2037
PTV liner, sintered glass, deactivated	1.5	112	5190-1426
Liners for High Temperature PTV Inlet, G3506A			
PTV liner, high temperature, quartz	3.4	713	5188-5313
PTV liner, high temperature, borosilicate	3.4	668	5188-5356

Syringes for Septumless and High Temperature PTV Inlets

Volume (μL)	Description	Needle	Part No.
0.5	Removable	23/70/HP	5182-9651
5	Straight, fixed	23/42/HP	9301-0892
10	Straight, fixed	23/42/HP	9301-0713
50	Straight, fixed, for large volume injections	23/42/HP	5183-0318
100	Straight, fixed, for large volume injections	23/42/HP	5183-2058

Purged Packed Inlets

Packed column analysis is frequently done when high efficiency separations are not needed or when gases are analyzed by gas-solid chromatography. Purged packed inlets are simple in both design and use. Few parameters need to be set, and all carrier gas flow flushes through the inlet into the column in the standard configuration.




7890/6890/6850 Purged Packed Inlet Supplies

Description	Unit	Part No.
QuickPick Purged Packed Inlet PM Kit		5188-6498
Includes 5 non-stick BTO septa, 1 O-ring, 1 ferrule, and 1 disposable glass liner		
Merlin Microseal		5182-3444
Microseal high pressure nut		5182-3445
Septum retainer nut		18740-60835
11 mm Certified BTO septa	50/pk	5183-4757
Top insert weldment		19243-80570
O-ring, Fluorocarbon	12/pk	5080-8898
Disposable glass liner, 170 μ L internal volume	25/pk	5080-8732
Disposable glass insert, deactivated, 170 μ L internal volume	5/pk	5181-3382
Ferrule, 1/4 in. Vespel	10/pk	5080-8774
1/4 in. nut, brass	10/pk	5180-4105
530 μ m column adapter for use with glass liners		19244-80540
1/8 in. column adapter for use with glass liners		19243-80530
1/4 in. column adapter for use with glass liners		19243-80540
Insulating cup		19234-60720
Universal column nut	2/pk	5181-8830

Nuts and Ferrules for 1/8 in. Packed Columns

Description	Unit	Part No.
1/8 in. stainless steel nut and ferrule set	20/pk	5080-8751
1/8 in. brass nut and ferrule set	20/pk	5080-8750
Vespel/graphite ferrule, 1/8 in.	10/pk	0100-1332



Flame Ionization Detector (FID)


FID Jet Identification and Selection

Before ordering parts for FID maintenance, determine which type of FID is installed on your GC. The FID is available in two versions:

- Dedicated, Capillary Optimized: for capillary columns only
- Adaptable: for packed or capillary columns


To determine the type of FID installed on your GC, open the oven door and examine the fitting at the base of the detector. Compare to the following diagram.

Hint: Adaptable jets are longer than dedicated capillary jets.



FID Jets

Description	Jet Tip ID	Length (mm)	Part No.
Jets for capillary optimized fittings			
Capillary	0.29 mm (0.011 in.)	42.8	G1531-80560
Capillary, high temperature Use with simulated distillation	0.47 mm (0.018 in.)	42.8	G1531-80620
Jets for adaptable fittings			
Capillary	0.29 mm (0.011 in.)	61.5	19244-80560
Capillary, high temperature Use with simulated distillation	0.47 mm (0.018 in.)	61.5	19244-80620
Packed	0.46 mm (0.018 in.)	63.5	18710-20119
Packed, wide-bore Use with high-bleed applications	0.76 mm (0.030 in.)	63.5	18789-80070



Flame Ionization Detector (FID) assembly

7890/6890/6850 Flame Ionization Detector (FID) Supplies

Item	Description	Unit	Part No.
	PTFE chimney (optional)		19231-21050
1	Collector nut		19231-20940
2	Spring washer	10/pk	5181-3311
3	Ignitor castle		19231-20910
	Hastelloy ignitor castle (optional)		19231-21060
4	Ignitor glow plug assembly		19231-60680
5	Collector insulator		G1531-20700
6	Collector body		G1531-20690
	Hastelloy collector body		G1531-21090
7	FID collector assembly		G1531-60690
	FID collector cleaning brush	2/pk	8710-1346
	Collector housing		G1531-20740
	FID retainer nut wrench 5880, 5890, 6890		19301-00150
	1/4 in. nut driver for FID jet, drilled shaft		8710-1561
8	FID ignitor cable for 6890/6850 only		G1531-60680
	FID ignitor cable, 7890A only		G3431-60680
	FID performance evaluation sample kit This sample is used for the HP 5880, 5890 and 6890 with a FID or TCD. Solution of 0.033% C14, C15, and C16 normal alkanes in hexane. Three 0.5 mL ampoules.		18710-60170
	FID MDL test sample for 7890 only 3 x 0.5 mL ampoules. Contains 2.36 mg/L n-Tridecane, 2.36 mg/L n-Tetradecane, 23.6 mg/mL n-Penta-decane, 23.6 mg/mL n-Hexadecane in iso-octane		5188-5372

(Continued)




FID cleaning kit, 9301-0985

7890/6890/6850 Flame Ionization Detector (FID) Supplies

Item	Description	Unit	Part No.
	O-rings	12/pk	5080-4978
	FID/NPD adapter for capillary column		19244-80610
	FID/NPD 1/8 in. packed column		19231-80520
	FID/NPD 1/4 in. packed column		19231-80530
	1/4 in. nut driver for FID jet, drilled shaft		8710-1561
	FID collector cleaning brush	2/pk	8710-1346
	FID Supplies Kit, Includes:		5182-3450
	Jet, packed standard 0.018 in. ID tip	3 each	18710-20119
	FID and TCD sample	2 each	18710-60170
	Ignitor glow plug assembly	2 each	19231-60680
	Jet, 0.011 in. ID tip, capillary adaptable	3 each	19244-80560
	FID flow measuring insert	2 each	19301-60660
	Cleaning wires for 0.03 in. ID jet	5/pk	5180-4150
	Cleaning wire for 0.018 in. ID/530 µm jet	5/pk	5180-4152
	For use with 0.018 and 0.011 in. ID jets		
	Cleaning wire for 0.011 in. ID jet	5/pk	19301-20720
	Capillary inlet cleaning wires	5/pk	5180-4153
	FID cleaning kit		9301-0985

Electron Capture Detector (ECD)

The Agilent micro ECD is the most sensitive on the market, with a detection zone volume 10 times smaller than any other ECD. The replaceable liner serves as a physical stop for the column, ensuring reproducible column installation and decreasing column contamination of the cell.



Electron Capture Detector (ECD) assembly

7890/6890/6850 Electron Capture Detector (ECD) Supplies

Item	Description	Part No.
1	Standard ECD makeup gas adapter*	G1533-80565
	Micro ECD makeup gas adapter	G2397-80520
	Micro ECD mixing liner, also compatible with standard ECD design	G2397-20540
	Gigabore liner for standard ECD, polyamide coating, not compatible with micro ECD	19233-20625
	ECD makeup gas adapter, 7890 only	G3433-63000
2, 3	Insulating cup	19234-60720
	ECD adapter end cap	19233-20755
	Vespel ferrule, 1/4 in.	5080-8774
	1/4 in. nut, brass	5180-4105
	Electron Capture Detector sample	18713-60040
	ECD test sample	5183-0379
	Micro ECD wipe test kit	18713-60050


*Includes one each of P/N 19233-20625 and 19233-20755

Thermal Conductivity Detector (TCD)


The TCD compares the thermal conductivities of two gas flows — pure carrier gas (also called the reference gas) and carrier gas plus sample components (also called column effluent).




1/8 in. SS packed column



1/4 in. SS packed column



Standard design



7890/6890/6850 Thermal Conductivity Detector (TCD) Supplies

Item	Description	Unit	Part No.
For 1/8 in. SS Packed Column Installation			
1	Vespel/graphite ferrule, 1/8 in.	10/pk	0100-1332
2	1/8 in. nut, brass	10/pk	5180-4103
For 1/4 in. SS Packed Column Installation			
3	Vespel/graphite ferrule, 1/8 in.	10/pk	0100-1332
4	1/8 in. nut, brass	10/pk	5180-4103
5	1/4 in. packed column adapter		G1532-20710
6	Vespel ferrule, 1/4 in.	10/pk	5080-8774
7	1/4 in. nut, brass	10/pk	5180-4105
For Capillary Column Installation (Standard)			
	TCD capillary column adapter		G1532-80540
8	Vespel/graphite ferrule, 1/8 in.	10/pk	0100-1332
9	1/8 in. nut, brass	10/pk	5180-4103
	Universal column nut	2/pk	5181-8830
	6850 column nut	2/pk	5183-4732
	530 μ m, 1.0 mm ID graphite ferrule	10/pk	5080-8773
	320 μ m, 0.5 mm ID graphite ferrule	10/pk	5080-8853
	TCD sample		18711-60060
	FID and TCD sample		18710-60170

TCD Ferrules

Column ID (mm)	Back Ferrules	Front Ferrules 10/pk
0.53	5182-3477	5182-9673
0.32	5182-3477	5182-9676
0.25/0.2/0.1	5182-3477	5182-9677
No hole	5182-3477	5182-9679
TCD back ferrule for 1/8 in. detector fitting, 10/pk		5180-4103

Determining the TCD Electronic Pressure Control (EPC)

If you have a 6890A or 6890A Plus GC, you may have an older design EPC flow manifold for the TCD. The older design requires removal of sheet metal panels to attach the TCD reference flow gas supply inside the GC. The new "Minifold" design allows TCD reference gas to be connected directly to the back of the GC. Replacement TCD filament block assemblies have different part numbers depending on the EPC design type.

Once you have determined the type of EPC module, consider ordering a passivated filament block assembly, which is recommended for fatty acid analysis or reactive/acidic samples.

TCD Filament Block Assemblies


Instrument	Passivated	Applications	Specifications	EPC Design	Part No.
7890A	Yes	Standard TCD Analysis Gases/Hydrocarbons	Complete Detector Assembly Includes detector palette and heater/sensor assembly	Original	G3432-60220
7890A	Yes	Standard TCD Analysis Gases/Hydrocarbons	Complete Detector Assembly Includes detector palette and heater/sensor assembly Third detector, side mounted	Original	G3432-60221
6890	No	Standard TCD Analysis Gases/Hydrocarbons	Filament Block Only Must reuse heater/sensor	Original	G1532-60675
6890	No	Standard TCD Analysis Gases/Hydrocarbons	Filament Block Only Must reuse heater/sensor	Minifold	G1532-60685
6890	Yes	Recommended for Fatty Acid Analysis	Filament Block Only Must reuse heater/sensor	Original	G1532-60690
6890/6850	Yes	Recommended for Fatty Acid Analysis	Filament Block Only Must reuse heater/sensor	Minifold	G1532-60695
6890/6850	No		Complete Detector Assembly Includes detector palette and heater/sensor assembly	Minifold	G2630-61230

Flame Photometric Detector (FPD)

In 2005, Agilent released an improved FPD with minimum detectable levels (MDL) of 3.6 pg/s for sulfur and 60 fg/s for phosphorus. This is more than a 5 times improvement for sulfur. The updated design is based on a one-piece, deactivated transferline jet assembly and improved optics. Upgrade kits are available.

7890/6890/6850 FPD Ignitor and Heat Shield Assembly


Item	Description	Part No.
1	FPD Exit Tube Assembly	
	Aluminum	19256-60700
	Stainless steel	19256-20705
2	Vespel ferrule, 1/4 in., 10/pk	5080-8774
3	Emission chamber	
	FPD, single	19256-80560
	FPD, dual	19256-80600
5	FPD ignitor replacement kit	19256-60800
	O-ring, size 2-010	0905-1610
	Spacer, ignitor	19256-20590
	Glow plug	0854-0141
5	Ignitor cable assembly	G1535-60600
6	Heat shield gasket, white	19256-80045
7	First heat shield window	19256-80030
8	Heat shield disk	19256-20580
9	Stainless steel coupling	19256-20550
10	Lock washer (4 required)	2190-0584
11	Screw, M3 x 12 mm, T10 (4 required)	0515-1084
12	Collar	19256-20690
13	Screw, M3 x 66 mm, T10	0515-0680
	FPD check out sample	5188-5953
	FPD sample	5188-5245



FPD ignitor and heat shield assembly

FPD Lens Assembly

Item	Description	Part No.
1	Clamp	19256-00090
2	Screw, M3 x 25 mm (4 required)	0515-0683
3	Window O-ring, inner, 0.926 in. ID, orange	5061-5886
4	Convex lens	1000-1438
5	Lens housing	19256-20900
6	Flange ring	19256-00200
7	Fluorocarbon Elastomer O-ring, brown, 1.239 in. ID	5061-5890



FPD lens assembly

FPD Photomultiplier Tube (PMT) and Bracket Assemblies

Description	Part No.
Chimney back cover	G1535-80520
Heator/sensor assembly	G1535-60610
Transfer line support bracket	19256-00320
Bracket/support	G1535-00010
Sulfur filter, 7890 and late model 6890*	1000-1437
Sulfur filter, blue, early model 6890*	19256-80000
Phosphorus filter, yellow	19256-80010
Filter spacer (used only with sulfur filter)	19256-20910
PMT housing assembly	19256-60510
Dual FPD chimney front	G1535-00030

*Please contact Agilent technical support for assistance in selecting the correct sulfur filter for your 6890 FPD detector.



BloS NPD bead assembly, G3434-60806

Nitrogen Phosphorus Detector (NPD)

NPD Beads

The NPD for the 7890/6890 GC features a ceramic bead selective for nitrogen and phosphorous compounds. Agilent offers three beads:

- BloS bead
- White ceramic bead
- Black ceramic bead

Compared to the white ceramic bead, the new BloS bead provides:

- Superior bead lifetime
- Faster attainment of stable operation at initial start-up, as well as more stable operation throughout bead's lifetime
- Superior sensitivity and selectivity for phosphorous-containing compounds
- Similar sensitivity and selectivity for nitrogen-containing compounds
- Superior immunity to moisture

The white ceramic bead exhibits some tailing for phosphorous compounds. The black ceramic bead does not exhibit peak tailing and typically has a longer lifetime than the white bead; however, it is less sensitive.

All Agilent NPD beads are preconditioned, self-aligning for installation and include a proof-of-performance chromatogram.

NPD Beads


Description	Part No.
BloS NPD bead assembly	G3434-60806
NPD white ceramic bead assembly	G1534-60570
NPD black ceramic bead assembly	5183-2007

NPD Jet Identification and Selection

Before ordering parts for NPD maintenance, determine which type of NPD is installed on your GC. The NPD is available in two versions:


- Dedicated, Capillary Optimized: for capillary columns only
- Adaptable: for packed or capillary columns

Hint: Adaptable jets are longer than dedicated capillary jets.



NPD Jets

Description	Jet Tip ID	Length (mm)	Part No.
Jets for capillary optimized fittings			
Capillary with extended jet (recommended)	0.29 mm (0.011 in.)	51.5	G1534-80580
Capillary	0.29 mm (0.011 in.)	42.8	G1531-80560
Capillary, high temperature	0.47 mm (0.018 in.)	42.8	G1531-80620
Jets for adaptable fittings			
Capillary with extended jet (recommended)	0.29 mm (0.011 in.)	70.5	G1534-80590
Capillary	0.29 mm (0.011 in.)	61.5	19244-80560
Capillary, high temperature	0.47 mm (0.018 in.)	61.5	19244-80620
Packed	0.46 mm (0.018 in.)	63.5	18710-20119



Nitrogen Phosphorus Detector (NPD) assembly

7890/6890/6850 Nitrogen Phosphorus Detector (NPD) Supplies

Item	Description	Part No.
1	Screws, M3 x 0.5, 8 mm (Pozidriv)	0515-0655
2	NPD white ceramic bead assembly NPD black ceramic bead assembly	G1534-60570 5183-2007
3	Lid weldment	G1534-80510
4	Metal C-ring, top	0905-2580
5	Alumina insulator, upper	G1534-40020
6	Collector funnel	G1534-20530
7	Alumina insulator, lower	G1534-40030
8	Metal C-ring, bottom	0905-1284
9	Screw, M4 x 07, 10 mm	0515-2495
10	J-Clamp	1400-0015
11	Screw, M4 x 07, 10 mm	0515-2495
12	NPD interconnect assembly	G1534-60610
13	Mounting pallet	G1531-40020
14	Jet, 0.011 in./0.29 mm ID tip, capillary dedicated Jet, 0.011 in. ID tip, capillary adaptable Jet, packed standard 0.018 in. ID tip	G1531-80560 19244-80560 18710-20119
15	Base weldment, capillary NPD for 6890/6850 only Base weldment, packed NPD for 6890/6850 only Base weldment, capillary NPD, 7890A Base weldment, packed NPD, 7890A	G1534-80500 G1534-80540 G3434-67500 G3434-67540
16	Lid stop NPD ceramic insulator kit Includes items 4, 5, 7, and 8	G1534-20590 5182-9722
17	FID/NPD adapter for capillary column FID/NPD 1/8 in. packed column FID/NPD 1/4 in. packed column	19244-80610 19231-80520 19231-80530
18	Insulating cup Vespel ferrule, 1/4 in., 10/pk 530 µm, 1.0 mm ID graphite ferrule, 10/pk 320 µm, 0.5 mm ID graphite ferrule, 10/pk 1/4 in. nut, brass, 10/pk Universal column nut, 2/pk Nitrogen Phosphorus Detector sample	19234-60720 5080-8774 5080-8773 5080-8853 5180-4105 5181-8830 18789-60060



Nitrogen Chemiluminescence Detector (NCD)

Nitrogen and Sulfur Chemiluminescence Detectors

The Agilent 355 Sulfur Chemiluminescence Detector (SCD) is the most sensitive and selective chromatographic sulfur detector available for the analysis of sulfur compounds.

The Agilent 255 Nitrogen Chemiluminescence Detector (NCD) is a nitrogen-specific detector that produces a linear and equimolar response to nitrogen compounds based on a chemiluminescent reaction of NO with ozone. Even complex sample matrices can be analyzed with little or no interference.



Dual plasma burner accessory kit, G6600-60038



Replacement oil coalescing filter, G6600-80042



Replacement oil coalescing filter for oil mist filter, G6600-80044



Replacement odor filtration element, G6600-80045

Nitrogen Chemiluminescence Detector (NCD) Supplies

Description	Part No.
Dual plasma burner accessory kit Includes ferrules, fittings and quartz tube	G6600-60038
PM Kit, dry piston pump Includes dry piston seal, 6 Moleculite vacuum pump traps and 1 set of ceramic tubes	G6600-67006
PM Kit, DP RV5 oil pump Includes 6 chemical traps for ozone destruction, 3 oil coalescer elements and 4 (1 qt) bottles of synthetic oil	G6600-67007
PM Kit, dry piston pump Includes 4 chemical traps for ozone destruction and 2 repair kits for pump	G6600-67008
12-month maintenance kit Includes 6 qt oil, 12 chemical traps, 4 oil filter elements, 2 sets of ceramic tubes, and 2 O-rings	G6600-67009
Replacement oil coalescing filter	G6600-80042
Oil mist filter for RV5 pump	G6600-80043
Replacement oil coalescing filter for oil mist filter	G6600-80044
Replacement odor filtration element	G6600-80045
O-ring, 1.3614 in. ID	G6600-80050
O-ring, 1.301 in. ID	G6600-80051
Dual plasma quartz tube	G6600-80063
Mobil 1 synthetic oil	G6600-85001
Oil, Edwards Ultragrade for RV3 and RV5 pumps	G6600-85002
Spare column nut and ferrule kit	G6600-80018



Sulfur Chemiluminescence Detector (SCD)



PM kit, G6600-67008



Dual plasma burner accessory kit, G6600-60037



Oil mist filter, G6600-80043

Sulfur Chemiluminescence Detector (SCD) Supplies

Description	Part No.
PM Kit, dry piston pump	G6600-67006
Includes dry piston seal, 6 Moleculite vacuum pump traps and 1 set of ceramic tubes	
PM Kit, dry piston pump	G6600-67008
Includes 4 chemical traps for ozone destruction and 2 repair kits for pump	
12-month maintenance kit	G6600-67009
Includes 6 qt oil, 12 chemical traps, 4 oil filter elements, 2 sets of ceramic tubes, and 2 O-rings	
Dual plasma burner accessory kit	G6600-60037
Mobil 1 synthetic oil	G6600-85001
Oil mist filter for RV5 pump	G6600-80043
Oil, Edwards Ultragrade for RV3 and RV5 pumps	G6600-85002
O-ring, 1.301 in. ID	G6600-80051
Ozone destruction chemical trap	G6600-85000
Replacement oil coalescing filter for oil mist filter	G6600-80044
Sulfur chemiluminescence test sample	G2933-85001
Sulfur trap	G2933-85003
For carrier H ₂ and air gases; one required for each cylinder of gas (3 total)	
Spare column nut and ferrule kit	G6600-80018

Miscellaneous Instrument Parts and Supplies

Description	Part No.
Oven exhaust deflector for 6890	G1530-80650
Oven exhaust deflector for 6850	G2630-60710

GC Standards

GC Qualitative Standards

Description	Part No.
Qualitative Simulated Distillation Standards	
Boiling Point Calibration Sample No. 1	5080-8716
Low Boiling Point Calibration Sample No. 220	5080-8768
Boiling Point Calibration Sample No. 320	5080-8769
PolyWax 500, 1 g, neat	5188-5316
PolyWax 655, 1 g, neat	5188-5317
Qualitative Petrochemical Standards	
Alcohol in Gasoline Sample	18900-60640
Natural Gas Sample	5080-8756
Transformer Gas Sample	5080-8759
Refinery Gas Sample	5080-8755
Reference Gas Oil No. 1, Batch 2	5060-9086
Miscellaneous Qualitative Standards	
Nickel Catalyst Test Sample	19354-60510
Nickel Catalyst refill	5080-8761
MIDI System Calibration Standard	19298-60500
Oral Fluids Analyzer Test Sample	G1540-85010



GC/MS Parts and Supplies

Your mass spectrometer is a sensitive, specialized device that delivers a higher level of functionality than other GC detectors. To continue achieving optimal results, it is critical to maintain your system properly by performing the essential tasks within this section. Some of the benefits of maintaining your GC/MSD include:

- Less downtime for repairs
- Longer lifetime for your MSD system
- Reduction in overall operating costs

It is advisable to keep a log book of system performance, Autotune, and maintenance operations performed. This makes it easier to identify variations from normal performance and to take corrective action.

Maintenance Schedule

Task	Every week	Every 6 months	Every year	As needed
Tune the MSD				◆
Change injection port liners	◆			
Check the foreline pump oil level	◆			
Gas ballast the foreline pump	◆			
Check the calibration vial		◆		
Replace the foreline pump oil		◆		
Check the diffusion pump fluid	◆			
Replace the diffusion pump fluid			◆	
Replace the dry pump diaphragm seals (MVP55)				◆
Replace the dry pump tip seals (IDP3)			◆	
Replace the traps and filters			◆	
Clean the ion source				◆
Change the carrier gas trap(s) and purifier				◆
Replace worn out parts				◆
Lubricate seals (where appropriate)				◆
Replace column				◆



For in-depth information about maintaining your GC/MS, request "*Maintaining Your Agilent GC and GC/MS Systems*" from your Agilent Representative (**publication number 5990-5451EN**).



MSD Contamination


Contamination is usually identified by excessive background in the mass spectra, which can come from the GC or MSD. The source of contamination can sometimes be determined by identifying the contaminants. Some contaminants are much more likely to originate in the GC, while others are likely to originate in the MSD.

MSD Contamination Identification

The following table lists some of the more common contaminants, the ion characteristics of those contaminants, and the likely sources of those contaminants.

Tips & Tools

A crude sign of a "leak-free" MS system is when the ion ratio of m/z 28 (nitrogen) over m/z 32 (oxygen) is approximately two or greater.



Common Contaminants

Ions (m/z)	Compound	Possible Source
13, 14, 15, 16	Methane	CI gas
18, 28, 32, 44 or 14, 16	H ₂ O, N ₂ , O ₂ , CO ₂ , CO or N, O	Residual air and water, air leaks, outgassing from Vespel ferrules
31, 51, 69, 100, 119, 131, 169, 181, 214, 219, 264, 376, 414, 426, 464, 502, 576, 614	PFTBA and related ions	PFTBA (tuning compound)
31	Methanol	Cleaning solvent
43, 58	Acetone	Cleaning solvent
78	Benzene	Cleaning solvent
91, 92	Toluene or Xylene	Cleaning solvent
105, 106	Xylene	Cleaning solvent
151, 153	Trichloroethane	Cleaning solvent
69	Foreline pump fluid or PFTBA	Foreline pump oil vapor or calibration valve leak
73, 147, 207, 221, 281, 295, 355, 429	Dimethylpolysiloxane	Septum bleed or methyl silicone column coating
77, 94, 115, 141, 168, 170, 262, 354, 446	Diffusion pump fluid	Diffusion pump fluid and related ions
149	Plasticizer (phthalates)	Vacuum seals (O-rings) damaged by high temperatures, use of vinyl or plastic gloves
Peaks spaced 14 amu apart	Hydrocarbons	Fingerprints, foreline pump oil

Cleaning and Maintenance Supplies

Description	Part No.
One Year Maintenance Kit (for diffusion pump systems) Includes Big Universal Trap for He (1/8 in.), abrasive sheets (5/pk), lint-free cloths (15/pk), cotton swabs (100/pk), SantoVac Ultra, 18.5 mL (2 each), rough pump oil (1 L), filament assembly, octafluoronaphthalene (OFN)	5183-2096
Nylon gloves, lint-free, large, 1 pair	8650-0030
Nylon gloves, lint-free, small, 1 pair	8650-0029
Lint-free industrial wipes, 100% cotton, 9 x 9 in., 300/pk	9310-4828
Ion source cleaning kit Includes lint-free cloths (15/pk), abrasive sheets (5/pk), cotton swabs (100/pk), lint-free nylon gloves, abrasive Alumina powder	5181-8863
Cloths, lint-free, 15/pk	05980-60051
Cotton swabs, 100/pk	5080-5400
Abrasive sheets, aluminum oxide green lapping paper, 600 mesh, 5/pk	5061-5896
Alumina powder, abrasive, 1 kg	8660-0791
PFTBA sample, certified, 10 g, 5.32 mL	8500-0656
Replacement glass bulb for PFTBA and PFDTD test sample	G3170-80002
Replacement glass vial for PFTBA and PFDTD test sample	05980-20018
Activated alumina, absorbent pellets for Edwards rough pump traps, non-LC/MS, 1 lb can	8500-1233
MSD Tool Kit, 5975/5973 Includes source hold tool, lint-free cloth, cotton swabs, lint-free nylon gloves, abrasive sheets, wrenches and driving tools	G1099-60566
MSD Tool Kit, 5972/5971 Includes small cleaning rod, large cleaning rod, source hold tool, cotton swabs, lint-free nylon gloves, abrasive sheets, wrenches and driving tools	05971-60561
MS Interface Supplies	
MS interface column nut, female	05988-20066
Inlet column nut for long or long two-hole ferrules	05921-21170
Universal column nut, 2/pk	5181-8830
MS interface column installation tool for 5973 and 5975	G1099-20030
Column installation tool for 5975T	G3880-20030

(Continued)



MS interface column nut, 05988-20066



Universal column nut, 5181-8830



Column installation tool, G1099-20030




Vespe/Graphite ferrules, 5181-3323

Cleaning and Maintenance Supplies

Description	Part No.
Tools	
Screwdriver, 3 in. Pozidriv shaft No. 1 pt, fits no. 2-4 screws	8710-0899
Screwdriver, 4 in. Pozidriv shaft No. 2 pt, fits no. 5-10 screws	8710-0900
Open end wrench, 1/4 and 5/16 in.	8710-0510
Hex nut driver, 5.5 mm	8710-1220
Screwdriver, Torx T20	8710-1615
Screwdriver, Torx T15	8710-1622
Screwdriver, Torx T10	5182-3466
Ferrules	
0.4 mm Vespe/Graphite ferrule for 200/250 μ m columns, 10/pk	5062-3508
0.5 mm Vespe/Graphite ferrule for 320 μ m columns, 10/pk	5062-3506
250 μ m Vespe/Graphite ferrule, 10/pk	5181-3323
SilTite metal ferrules for 1/16 in. OD tubing, 10/pk Includes 2 column nuts	5184-3571
SilTite metal ferrules, 1/16 in. x 0.4 mm ID, 10/pk Includes 2 column nuts	5184-3569
SilTite metal ferrules, 1/16 in. x 0.5 mm ID, 10/pk Includes 2 column nuts	5184-3570
Ferrule pre-swaging tool	G2855-60200
Plug for microfluidic manifold or unions	G2855-60570

Tips & Tools

Even preconditioned ferrules can shrink slightly at very high temperatures. If leak problems persist upon a new column installation, check this fitting first.





Electron Impact (EI) Ion Source

Ion Source

The ion source operates by electron ionization (EI) or chemical ionization (CI). The sample enters the ion source from the GC/MSD interface. Electrons emitted by a filament enter the ionization chamber, guided by a magnetic field. The high-energy electrons interact with the sample molecules, ionizing and fragmenting them. The positive voltage on the repeller pushes the positive ions into the lens stack, where they pass through several electrostatic lenses. These lenses concentrate the ions into a tight beam, which is directed into the mass filter.


Electron Impact (EI) Ion Source

The recommended cleaning material for the EI ion source is abrasive, aluminum oxide powder.

Do not immerse filaments or lens insulators in solvent. If insulators are dirty, clean them with a cotton swab dampened with reagent-grade methanol. If that does not clean the insulators, replace them.

5975/5973 MSD Electron Impact Ion Source Parts (EI)


Item	Description	Part No.	Inert Part No.
1	Ion source body	G1099-20130	G2589-20043
2	Drawout plate, 3 mm	05971-20134	G2589-20100
	Drawout plate, 6 mm	G3163-20530	G2589-20045
3	Drawout cylinder	G1072-20008	G1072-20008
4	Ion focus lens	05971-20143	05971-20143
5	Lens insulator	G3170-20530	G3170-20530
6	Entrance lens	G3170-20126	G3170-20126
7	Cap screw, gold plated	G1999-20021	G1999-20021
8	High temperature filament	G2590-60053	G2590-60053
9	Transfer line socket	G1099-20136	G1099-20136
10	Repeller insulator	G1099-20133	G1099-20133
11	Repeller	G1099-20132	G2589-20044



5975/5973 MSD Electron Impact (EI) ion source assembly

Tips & Tools

It is good practice to replace scratched lenses and other ion source parts regularly. Scratched source parts lead to poor performance.

**5972/5971/GCD MSD Ion Source Parts (EI)**


Description	Part No.
Entrance lens	05971-20126
Lens insulator	G3170-20530
Ion focus lens	05971-20143
Drawout cylinder	G1072-20008
Drawout plate, 3 mm	05971-20134
Set screw	0515-1446
Repeller assembly	05971-60170
Screw for filament on the source	0515-1046
Transfer line tip, gold plated	05971-20305

Warnings & Caution

Important: Do not abrasively or ultrasonically clean the insulators.


Abrasively clean the surfaces that contact the sample or ion beam. Use an abrasive slurry of alumina powder and reagent-grade methanol on a cotton swab. Use enough force to remove all discoloration. Polishing the parts is not necessary; small scratches will not harm performance. Abrasively clean discoloration where electrons from filaments enter the source body.

Take care to avoid contaminating cleaned and dried parts. Put on new, clean gloves before handling the parts. Do not put the cleaned parts on a dirty surface. Place them only on clean, lint-free cloths.



Tips & Tools

Visual appearance is not an accurate guide to cleanliness of the CI ion source. The CI ion source can show little or no discoloration, yet still need cleaning.



Chemical Ionization (CI) Ion Source


Because the CI ion source operates at much higher pressures than the EI ion source, it will probably require more frequent cleaning than the EI ion source.

The source should be cleaned whenever there are performance anomalies that are associated with a dirty ion source. Let analytical performance be your guide.

When cleaning the CI ion source, concentrate on the CI repeller, ion source body, and drawout plate. Be sure to clean the 0.5 mm diameter holes in the ion source body and drawout plate.


Cleaning the ion source is very similar to cleaning the EI ion source. Use the same EI cleaning procedure with the following exceptions:

- The CI ion source may not look dirty, but deposits left by chemical ionization are very difficult to remove. Clean the CI ion source thoroughly.
- Use a round wooden toothpick to gently clean out the electron entrance hole in the source body and the ion exit hole in the drawout plate.
- Do not use halogenated solvents. Use hexane for the final rinse.



5975/5973 MSD Chemical Ionization Ion Source Parts (CI)

Item	Description	Part No.
1	Source body	G1999-20430
2	Drawout plate	G1999-20446
3	Drawout cylinder	G1999-20444
5	Ion focus lens	G1999-20443
4	Lens insulator	G3170-20540
6	Entrance lens	G3170-20126
7	Repeller insulator	G1999-20433
8	Repeller	G1999-20432
9	High temperature filament	G1099-80053
10	Cap screw, gold plated	G1099-20021
11	Interface tip seal/spring	G1999-60412



5975/5973 MSD Chemical Ionization (CI) ion source assembly

QuickSwap MS Interface Restrictors

Agilent's QuickSwap Capillary Flow Technology module and pre-swaged fused silica tubing restrictors can increase the productivity of your Agilent 5973N and 5975 Inert MSD Systems, allowing you to change columns without venting the MSD. QuickSwap not included.

These restrictors are prefabricated for convenience and ease-of-use. For applications requiring other restrictor sizes, Agilent offers a wide variety of deactivated fused silica tubing, SilTite ferrules and swaging tools.

QuickSwap MS Interface Restrictors

Description	ID (mm)	Unit	Part No.
QuickSwap restrictor	0.092	4/pk	G3185-60361
QuickSwap restrictor	0.100	4/pk	G3185-60362
QuickSwap restrictor	0.110	4/pk	G3185-60363
QuickSwap restrictor	0.120	4/pk	G3185-60364
QuickSwap restrictor variety pack, 2 each of the above ID restrictors			G3185-60300

SilTite Metal Ferrules

Description	Unit	Part No.
For use with 0.25 mm ID capillary columns	10/pk	5188-5361
For use with 0.32 mm ID capillary columns	10/pk	5188-5362
For use with 1/16 in. OD stainless steel tubing Includes 2 column nuts	10/pk	5184-3571
For use with 0.53 mm ID capillary columns	10/pk	5188-5363



Filament assembly (EI), G3170-60050

MSD Filaments

Like the filaments in an incandescent light bulb, the ion source filaments will eventually burn out. Certain practices will reduce the chance of early failure.

- When setting up data acquisition parameters, set the solvent delay so that the analyzer will not turn on while the solvent peak is eluting
- When the software prompts 'Override solvent delay at the beginning of a run' always select 'No'
- Higher emission current will reduce filament life
- If you control your MSD from the Edit Parameters screen, always select 'MS Off' before changing any of the filament parameters

MSD Filaments

Description	5975 Series	5975T Series	5973 Series	5972 Series	5971 Series
Filament assembly (EI)	G3170-60050		G3170-60050	G3170-60050	05971-60140
Filament assembly (CI)	G1099-80053		G1099-80053		
Micro ion vacuum gauge	G3170-80001				
Triode gauge tube for measuring vacuum			0960-0897		
Ion gauge controller		G3880-80010			
Ion gauge tube		G3880-80011			

Tips & Tools

It is very useful to switch from one filament to the other every three months so that when a filament fails, you know the other will fail soon. This will allow you to change both filaments at the same time. Since the GC/MS system is already vented, it's a good idea to replace other supplies in the flowpath at the same time as the filaments.





Electron multiplier replacement horn

MSD Electron Multipliers and Replacement Horn

The lifetime of an electron multiplier is directly related to the current that flows through it and the extent of contamination or condensation that it experiences. Replace the electron multiplier or replacement horn when voltage is over 2500 V. To maximize electron multiplier life:

- Maintain the best possible vacuum, especially in the analyzer manifold
- Use extreme caution and be conservative with venting, pumpdown, and all vacuum system procedures to keep pump fluid background to a minimum
- After venting, allow four hours for pumpdown and thermal equilibration before scanning
- Actively look for background contamination and leaks and repair them immediately
- Don't tune excessively – PFTBA can result in higher background over an extended period of time
- Replace the electron multiplier if vacuum is poor or voltage is over 2500 V

MSD Electron Multipliers and Replacement Horn

Description	5975 Series	5973 Series	5972/5971 Series
Electron multiplier replacement horn	05971-80103	05971-80103	05971-80103
Use with electron multipliers with "straight" horns			
Triple axis detector assembly*	G3170-80100		
Triple axis electron multiplier	G3170-80103		
EM signal wire, low noise detector	G3170-80008		
High energy dynode		G1099-80001	
Electron multiplier			05971-80102

*Included on 5975 triple axis detector systems

Tips & Tools

The Agilent multipliers and horns listed are recommended for your MSD. Other manufacturers' products may be incompatible with Agilent instruments and can result in reduced sensitivity, lifetime, and noise problems.





Foreline Pump

Vacuum Systems and Pumps

Diffusion Pump

It is not necessary to change the diffusion pump fluid more than once a year, unless you observe symptoms that suggest a problem with the fluid. The MSD must be vented in order to check the diffusion pump fluid (except for the 5975/5973). Therefore, the best time to check the fluid is when the instrument is already vented for other maintenance.

Foreline Pump

The oil in the foreline or rough pump should be replaced on average once every six months, but can vary depending upon applications. If a foreline trap is present, the molecular sieves should also be replaced after an oil change.

Avoid contact with the pump oil. The residue from some samples may be toxic.
Dispense of used oil properly.

Pump Oils

Description	Part No.
Foreline pump oil, Inland 45, 1 L	6040-0834
High vacuum grease, 25 g	6040-0289
Diffusion pump fluid, 18.5 mL	6040-0809*
IDP Series tip seal kit for 5975T	IDP3TS

*2 required for 5975 and 5973 Series



7000A Triple Quadrapole GC/MS Parts and Supplies

Engineered from the ground up for ease-of-use and routine high performance operation, the 7000A Triple Quadrapole GC/MS delivers advanced high-speed GC/MS/MS quantitation for ultra-trace analysis of even the most complex samples. Combined with the Agilent 7890 GC, the result is an optimally robust GC/MS/MS system.

Gas Filters

Description	Part No.
Chemical Ionization Gas Purifier	G1999-80410
Big Universal Trap, 1/8 in. fittings, Helium (Ar/Me)	RMSH-2
Big universal trap, 1/8 in. fittings, Nitrogen	RMSN-2
Mounting clip	UMC-2

Foreline Pump Supplies

Description	Part No.
Diffusion pump fluid, 18.5 mL	6040-0809
Foreline pump oil, P3, 0.5 L	6040-0621
Rough pump inlet flange	0905-1463
Oil return kit	3162-1057
Pump oil drip pan	G1946-00034
Oil mist exhaust filter	G1099-80039
Oil mist filter for RV5 pump	G6600-80043



Low noise EM horn, G3170-80103



Filament assembly (EI), G3170-60050

Maintenance Supplies

Description	Part No.
Abrasive sheets	5061-5896
Alumina powder, abrasive, 1 kg	8660-0791
Cloths, lint-free	05980-60051
Lint-free industrial wipes, 100% cotton	9310-4828
Cotton swabs	5080-5400
Nylon gloves, lint-free, large	8650-0030
Nylon gloves, lint-free, small	8650-0029
High vacuum grease, 25 g	6040-0289
Electron multiplier replacement horn	05971-80103
Low noise EM horn	G3170-80103
Filament assembly, high temperature (EI)	G3170-60050
Filament assembly (CI)	G1099-80053
Micro ion vacuum gauge	G3170-80001
Replacement glass bulb for PFTBA and PFDTD test sample	G3170-80002



240-MS Ion Trap Parts and Supplies

The Agilent 240-MS Ion Trap delivers unparalleled capabilities for both research and routine applications. Advanced ionization, including positive and negative chemical ionization, improves selectivity and limits of detection. Enhanced scanning techniques ensure compound confirmation. The MS/MS and MSⁿ reduce matrix influences and provide more detailed structural information. The software comes with a full complement of productivity, reporting, and regulatory compliance tools.

- Accurate identification and quantification of trace analytes
- Unsurpassed sensitivity (200 femtogram OFN full scan)
- Choice of internal or external ionization configurations
- Powerful MS/MS and CI options
- Low maintenance and high reliability
- Intuitive software for increased productivity

240-MS Ion Trap Parts and Supplies

Description	Part No.
Manifold O-ring	393010924
Transfer line inner O-ring	393010920
Transfer line outer O-ring	393010918
Internal filaments (2 filaments on one disk)	392017401
Internal transfer line tip	393171201
External filament (single filament)	393161001
Electrode, end cap, Silchrom	393164493
Electrode set kit, Silchrom, DFC (inert) tested	9300003590
Includes 2 end cap electrodes, 1 RF electrode, cleaning instructions	
Electrode, RF, Silchrom	393167593
Spacer, RF, Silco-quartz	393053502
Electron multiplier	393175101
Transfer line assembly upgrade field kit	393101291
Contains a complete transfer line and side-mounted block for vacuum manifold	
EPA volatile kit for EPA methods 524.2 & 8260B	393082491
ChromatoProbe microvials, 100/pk	392567111

240-MS Ion Trap Parts and Supplies

Description	Part No.
GC/MS Standards	
Evaluation standard (Internal EI & CI) 2 pg/μL OFN, 5 pg/μL	393112601
Test standard for external EI (5 pg/μL OFN)	393112702
Benzophenone CI sensitivity standard 50 pg/mL	392030500
Test standard for external NCI (1 pg /μL DFB)	393113001
Tuning calibration compound PFTBA (FC-43)	392035300
GC/MS column test mix	392027300
Vacuum Supplies	
Oil mist exhaust filter, DS42	393847701
Oil mist eliminator	2735000500
Replacement cartridge for oil exhaust filter, 2/pk	2710100200
Foreline (roughing) pump oil, 1 L	8829951700
Premium foreline (roughing) pump oil, 1 L	8829953800
IDP-3 dry scroll pump tip seal maintenance kit	2710100400
IDP-3 dry scroll replacement module	2710100500





220-MS Parts and Supplies

The 220-MS is a high sensitivity, flexible gas chromatograph/mass spectrometer that delivers outstanding qualitative and quantitative data in a range of applications. This simple and robust system is easy to operate and maintain.

- Accurately identify and quantify trace analytes
- Take advantage of powerful CI and MS/MS upgrades for advanced applications
- Spend less time on maintenance and more time on analysis

220-MS Parts and Supplies

Description	Part No.
Electron multiplier assembly	393031501
Exit end cap electrode, chrome	393050292
Exit end cap electrode, SilChrom	393050293
Filament end cap electrode, chrome	393050392
Filament end cap electrode, SilChrom	393050393
RF ring electrode, chrome	393050492
RF ring electrode, SilChrom	393050493
Complete set of SilChrom electrodes and Silco-quartz spacers	393001991
Spacer, RF, quartz	393053501
Spacer, RF, Silco-quartz	393053502
Filament disk assembly with wire connectors	393060191
Filament disk assembly	392043700
User must solder on 3 wire connectors	
Thermocouple vacuum gauge	2722990700
Mass spectrometer expendable supplies kit for 2x0MS	393011391
Includes PFTBA calibration compound, cal-gas glass chamber, capillary injector nut, O-rings, cotton tipped applicators, end cap insulator, vacuum pump oil	
GC/MS Standards	
Hexachlorobenzene EI sensitivity standard 100 pg/mL	392027500
Benzophenone CI sensitivity standard 50 pg/mL	392030500
Tuning calibration compound PFTBA (FC-43)	392035300
Hexachlorobenzene EI sensitivity standard 2 pg/mL	392047100
GC/MS column test mix	392027300



MS standards

GC/MS Standards

GC/MS Analyzer Kit Standards


Description	Part No.
GC/MS semivolatiles analyzer checkout mixture	5190-0473
GC/MS pesticide analyzer internal standard	5190-0472
Pesticide analyzer checkout solution	5190-0468
Pesticide checkout standard, 100 µg/L, 3 x 1 mL	5190-0494
GC/MS toxicology checkout mixture	5190-0471
Residual solvent revised method 467, class C	5190-0493
Residual solvent revised method 467, class 1	5190-0490
Butanetriol internal standard #1 for biodiesel	5982-0024
Tricaprin internal standard #2 for biodiesel	5982-0025

MS Test and Performance Samples

	Description	Part No.	5975 Series	5973 Series	5972 Series	5971 Series	GCD	7000 Series
Tuning Samples								
El Tune	PFTBA sample, certified, 10 g, 5.32 mL	8500-0656	◆	◆	◆	◆	◆	◆
CI Tune	PFTBA MS Sample Kit, 0.942 g, 0.5 mL	05971-60571	◆	◆			◆	◆
	PFTD calibrant	8500-8510	◆	◆				◆
Performance Verification Samples								
EI	OFN, 1 pg/μL	5188-5348	◆	◆				
	Hexachlorobenzene 10 pg/μL, 1 ng/μL	8500-5808			◆			
	Methyl stearate (in methanol); 1 ng/μL, 2 ea	05990-60075				◆		
	Sample A, 10 ng/μL	05970-60045					◆	◆
Negative Mode CI	OFN, 100 fg/μL	5188-5347	◆					
Positive Mode CI	Benzophenone, 100 pg/μL	8500-5440	◆	◆	◆	◆		◆
Checkout Samples								
HighMass	PHFT, 100 pg/μL	5188-5357	◆					
Semi-Volatile	GC/MS tuning standard, DFTPP	8500-5995	◆	◆	◆	◆	◆	
Volatile	p-Bromofluorobenzene (BFB), 25 μg/mL	8500-5851	◆	◆	◆	◆	◆	
Evaluation sample	Solution of dodecane, biphenyl, p-chlorodiphenyl, and Methyl palmitate in isooctane. Six 1.0 mL ampoules: 4 at 10 ng/μL, 1 at 100 ng/μL, 1 at 100 pg/μL.	05970-60045	◆	◆	◆	◆		

Tips & Tools


Each GC/MS has a specific test and performance sample. Refer to the chart above for the exact sample. All volumes are approximately 0.5-1 mL unless otherwise specified.



GC Parts and Supplies for Varian Instruments


Agilent now provides replacement parts and supplies for GC instruments formerly manufactured by Varian. The following section includes ordering information for GC injectors, detectors and GC/MS systems.

Injector Replacement Parts and Supplies



1177 Split/Splitless Injector


Item	Description	Part No.
1	Injector nut	392597501
	Injector nut wrench	390842300
2	Knob	392597101
3	Automatic start switch	390820601
4	Septum, 9 mm	
	BTO	8010-0217
	Long Life	8010-0233
	Advanced Green	8010-0201
	Septum pick	7200008400
5	Septum purge head	
	EFC21 (stainless steel)	392597301
	EFC21 (UltiMetal)	392597303
	EFC25 or Manual Pneumatics	392597302
6	Purge head screw	391866308
7	Graphite liner O-ring, splitless, 6.5 mm	8004-0202
	Viton liner O-ring, 6.3 mm	8004-0201
8	Glass liner	8004-0165
9	Injector body	
	Stainless steel	392599401
	UltiMetal	392599411
	Manual	392599501
10	For replacement ferrules, see Agilent's new line of CrossLab Supplies	
11	Bottom nut	8004-0311

**1079 Large Volume Injector (LVI)**


Item	Description	Part No.
1a	Injector nut	394966601
1b	Injector nut	394966601
	Injector nut wrench	390842300
2	Automatic start switch	390820601
3	Septum, 11.5 mm	
	BTO	8010-0225
	Long Life	8010-0241
	Advanced Green	8010-0209
	Septum pick	7200008400
4	Septum support	391867600
5	Insert ferrule	8004-0204
6	Glass liner	8004-0164
7	Injector body, EFC type	
	Stainless steel	392544001
	UltiMetal	392544011
8	For replacement ferrules, see Agilent's new line of CrossLab Supplies	
9	Bottom nut	8004-0311

To learn more about Agilent CrossLab and to request your copy of the Agilent CrossLab product catalog, visit www.agilent.com/chem/CrossLab




**1093 Cool On-Column (COC) Injector**

Item	Description	Part No.
1	Injector nut	394966601
	Injector nut wrench	390842300
2	Automatic start switch	390820601
3	Septum, 11.5 mm	
	BTO	8010-0225
	Long Life	8010-0241
	Advanced Green	8010-0209
	Septum pick	7200008400
4	Septum support	391821100
5	Glass liner	
	Default	8004-0162
	High performance	8004-0167
	On-column	8004-0166
6	Screw	391866306
7	Graphite/Vespel ferrule	8004-0217
	Graphite ferrule	8010-0305
8	Bottom nut	
	Brass	8004-0311
	Stainless steel	8004-0312

**1061 Packed/530 µm Capillary Column Injector**

Item	Description	Part No.
1a	Injector nut	390812700
1b	Injector nut	392595501
	Injector nut wrench	390842300
2	Septum, 9.5 mm	
	BTO	8010-0219
	Long Life	8010-0235
	Advanced Green	8010-0203
	Septum pick	7200008400
3	Automatic start switch	390820601
4	Glass liner	8004-0168
5	Injector body, EFC23	392548301
6	Graphite/Vespel ferrule	8004-0217
	Graphite ferrule	8010-0305
7	For replacement ferrules, see Agilent's new line of CrossLab Supplies	
8	Bottom nut	8004-0311

**1041 Packed/Wide Bore On-Column (PWOC) Injector**

Item	Description	Part No.
1a	Injector nut	390812700
1b	Injector nut	392595501
	Injector nut wrench	390842300
2	Septum, 9.5 mm	
	BTO	8010-0219
	Long Life	8010-0235
	Advanced Green	8010-0203
	Septum pick	7200008400
3	Automatic start switch	390820601
4	Injector body, EFC type	392548201
5	Graphite/Vespel ferrule	8004-0217
	Graphite ferrule	8010-0305
6	Injector insert, stainless steel	392543101
7	For replacement ferrules, see Agilent's new line of CrossLab Supplies	
8	Bottom nut	8004-0311

Detector Replacement Parts and Supplies

Thermal Conductivity Detector (TCD)


Description	Part No.
Adapter TCD/DEFC capillary make-up	392585291
Adapter TCD/DEFC reference gas kit	392585292
Adapter TCD capillary make-up, MPC, 3800	392560591
TCD DEFC 14 (Non-H ₂), 2 channels	392561290

Flame Ionization Detector (FID)

Description	Part No.
Tube collector	394958700
Lower FID insulator 17311	2100003200
FID flame tip jet, 0.010 in.	200187500
FID flame tip jet with nut, 0.020 in.	200193800
Crunch washer, 25/pk	1500334701

Tips & Tools

To learn more about Agilent's complete portfolio of services, please visit www.agilent.com/chem/services



Pulsed Flame Photometric Detector (PFPD)

Description	Part No.
Photomultiplier tube (PFPD) #R647-08	392517100
O-Ring, silicone, 0.53 in. ID, PFPD	2740292400
PFPD light pipe	392515500
Sapphire window assembly	392514500
Sapphire window washer	392514300
Wrench, PFPD combustor support	392519200
Seal, combustor support	392513800
Combustor holder (2 mm)	392517800
Combustor Sulfur (2 mm), cleaned	392517600
Holder, combustor, 3 mm, cleaned	392517901
Combustor Phosphorus, 3 mm, cleaned	392517700

PFPD Filter Assemblies

Description	Part No.
Arsenic (As)	392515105
Manganese (Mn)	392544391
Nitrogen (N)	392511901
Sulfur and Phosphorus (S and P)	392515104
Phosphorus (P)	392515102
Sulfur (S)	392515101
Tin (Sn)	392515103

PFPD Nitrogen Mode Maintenance

Description	Part No.
Photomultiplier tube, Nitrogen R-5070A	392512800
O-Ring, 0.987 in. ID	2740236100
PFPD filter assembly, Nitrogen	392511901
PFPD light pipe	392515500
Sapphire window assembly	392514500
Sapphire window washer	392514300

Thermionic Specific Detector (TSD)

Description	Part No.
TSD bead probe, unconditioned and untested	390607400
TSD bead probe, conditioned and tested	390607401
Upper TSD insulator #17310 TSD	2100003100
O-Ring, 30/pk	2740928202
TSD collector assembly	390607900
Lower FID insulator #17311	2100003200
Crunch washer, 25/pk	1500334701
FID flame tip jet with nut, 0.020 in.	200193800
Flow tube assembly	200187600

GC/MS System Replacement Parts and Supplies

210/220/225 GC/MS Systems

Description	Part No.
Electron multiplier assembly	393031501
Exit end cap electrode, chrome	393050292
Exit end cap electrode, SilChrom	393050293
Filament end cap electrode, chrome	393050392
Filament end cap electrode, SilChrom	393050393
RF ring electrode, chrome	393050492
RF ring electrode, SilChrom	393050493
Complete set of SilChrom electrodes and Silco-quartz spacers	393001991
Spacer, RF, quartz	393053501
Spacer, RF, Silco-quartz	393053502
Filament disk assembly with wire connectors	393060191
Filament disk assembly	392043700
User must solder on 3 wire connectors	
Thermocouple vacuum gauge	2722990700
Mass spectrometer expendable supplies kit for 2x0MS	393011391
Includes PFTBA calibration compound, cal-gas glass chamber, capillary injector nut, O-rings, cotton tipped applicators, end cap insulator, vacuum pump oil	
GC/MS Standards	
Hexachlorobenzene EI sensitivity standard 100 pg/mL	392027500
Benzophenone CI sensitivity standard 50 pg/mL	392030500
Tuning calibration compound PFTBA (FC-43)	392035300
Hexachlorobenzene EI sensitivity standard 2 pg/mL	392047100
GC/MS column test mix	392027300

240 GC/MS and 4000 GC/MS Systems

Description	Part No.
Manifold O-ring	393010924
Transfer line inner O-ring	393010920
Transfer line outer O-ring	393010918
Internal filaments (2 filaments on one disk)	392017401
Internal transfer line tip	393171201
External filament (single filament)	393161001
Electrode, end cap, SilChrom	393164493
Electrode set kit, SilChrom, DFC (inert) tested	9300003590
Includes 2 end cap electrodes, 1 RF electrode, cleaning instructions	
Electrode, RF, SilChrom	393167593
Spacer, RF, Silco-quartz	393053502
Electron multiplier	393175101
Transfer line assembly upgrade field kit	393101291
Contains a complete transfer line and side-mounted block for vacuum manifold	
EPA volatile kit for EPA methods 524.2 & 8260B	393082491
ChromatoProbe microvials, 100/pk	392567111
GC/MS Standards	
Evaluation standard (Internal EI & CI) 2 pg/μL OFN, 5 pg/μL	393112601
Test standard for external EI (5 pg/μL OFN)	393112702
Benzophenone CI sensitivity standard 50 pg/mL	392030500
Test standard for external NCI (1 pg /μL DFB)	393113001
Tuning calibration compound PFTBA (FC-43)	392035300
GC/MS column test mix	392027300
Vacuum Supplies	
Oil mist exhaust filter, DS42	393847701
Oil mist eliminator	2735000500
Replacement cartridge for oil exhaust filter, 2/pk	2710100200
Foreline (roughing) pump oil, 1 L	8829951700
Premium foreline (roughing) pump oil, 1 L	8829953800
IDP-3 dry scroll pump tip seal maintenance kit	2710100400
IDP-3 dry scroll replacement module	2710100500

Saturn 2000 Series MS Systems

Description	Part No.
Mass spectrometer expendable supplies kit for 2x0MS Includes PFTBA calibration compound, cal-gas glass chamber, capillary injector nut, O-rings, cotton tipped applicators, end cap insulator, vacuum pump oil	393011391
Electron multiplier assembly	393031501
Filament disk assembly with wire connectors	393060191
Filament disk assembly User must solder on 3 wire connectors	392043700
Exit end cap electrode, SilChrom	393050293
Filament end cap electrode, SilChrom	393050393
RF ring electrode, SilChrom	393050493
Exit end cap electrode, chrome	393050292
Filament end cap electrode, chrome	393050392
RF ring electrode, chrome	393050492
Spacer, RF, quartz	393053501
Spacer, RF, Silco-quartz	393053502
Complete set of SilChrom electrodes and Silco-quartz spacers	393001991
GC/MS Standards	
Tuning calibration compound PFTBA (FC-43)	392035300
Benzophenone CI sensitivity standard 50 pg/mL	392030500
Hexachlorobenzene EI sensitivity standard 2 pg/mL	392047100
Hexachlorobenzene EI sensitivity standard 100 pg/mL	392027500
GC/MS column test mix	392027300

Tips & Tools

More information is a click away. We have a variety of educational primers, application notes, maintenance guides, and literature available from Agilent for free.

To learn more, visit www.agilent.com/chem/OnlineLibrary

